AutoSoft™ 30

Instructions for use / Mode d'emploi / Instrucciones de uso / Anwendungshinweise / Istruzioni per l'uso / Gebruiksaanwijzing / Brugsanvisning / Käyttöohje / Bruksanvisning / Instruktioner för Användning / Instrukcja użycia / Návod k použití / Οδηγίες χρήσης / Használati útmutató / הוראות השימוש / 使用说明 / إرشادات الاستعمال / Indicações de utilização / Инструкции по применению

TANDEM™
Diabetes Care

EN	Do not re-use
ES	No reutilizar
FR	Ne pas réutiliser
DE	Nicht wiederverwenden
IT	Non riutilizzare
NL	Niet hergebruiken
DA	Må ikke genbruges
FI	Ei saa käyttää uudelleen
NO	Må ikke brukes om igjen
SV	Får ej återanvändas
PL	Nie używać ponownie
CS	Nepoužívejte opakovaně
EL	Μην επαναχρησιμοποιείτε
HU	Nem újrafelhasználható
HE	לא לשימוש חוזר
ZH	请勿重复使用
AR	لا تُعيد استخدامه
PT	Não reutilizar
RU	Для одноразового использования

EN	Caution: see instructions for use
ES	Precaución: ver instrucciones de uso
FR	Attention : voir notice d'utilisation
DE	Achtung: siehe Gebrauchsanweisung
IT	Attenzione: consultare le istruzioni per l'uso
NL	Let op: zie gebruiksaanwijzing
DA	Forsigtig: Se brugsanvisningen
FI	Vaara: katso käyttöohjeet
NO	OBS! se bruksanvisningen
SV	Försiktighet: se bruksanvisningen
PL	Przeostoga: patrz instrukcja użycia
CS	Upozornění: viz návod k použití
EL	Προσοχή: ανατρέξτε στις οδηγίες χρήσης
HU	Vigyázat: olvassa el a használati útmutatót
HE	זהירות: יש לעיין בהוראות השימוש
ZH	注意事项: 见使用说明
AR	تحذير: راجع إرشادات الاستعمال
PT	Atenção: Ver as instruções de utilização
RU	Внимание: см. инструкцию по применению

EN	Consult instructions for use
ES	Consultar instrucciones de uso
FR	Consulter le mode d'emploi
DE	Anwendungshinweise beachten
IT	Consultare le istruzioni per l'uso
NL	Raadpleeg de gebruiksaanwijzing
DA	Se brugsanvisningen
FI	Lue käyttöohje
NO	Se bruksanvisningen
SV	Se bruksanvisningen
PL	Zapoznać się z instrukcją użycia
CS	Přečtete si návod k použití
EL	Συμβουλευτείτε τις οδηγίες χρήσης
HU	Olvassa el a használati útmutatót
HE	יש לעיין בהוראות השימוש
ZH	查阅使用说明
AR	راجع إلى إرشادات الاستعمال
PT	Consultar as instruções de utilização
RU	См. инструкцию по применению

EN	Prescription required, does not apply in all geographies
ES	La prescripción del producto no es aplicable en todas las áreas geográficas
FR	Soumis à prescription (ne concerne pas tous les marchés)
DE	Verschreibungspflichtig, gilt nicht für alle Regionen
IT	Prescrizione medica richiesta, non si applica in tutte le aree geografiche
NL	Uitsluitend op voorschrift, niet in alle landen van toepassing
DA	Receptpligtig; gælder ikke i alle lande
FI	Resepti vaaditaan; ei koske kaikkia alueita
NO	Resept påkrevd, gjelder ikke i alle geografiske områder
SV	Receptbelagt läkemedel, gäller inte i alla områden.
PL	Produkt dostępny na receptę; nie dotyczy niektórych lokalizacji
CS	Pouze na lékařský předpis. Neplatí ve všech zemích
EL	Απαιτείται συνταγή, δεν ισχύει σε όλες τις περιοχές
HU	Receptre kapható (nem minden országban)
HE	נדרש מרשם, לא חל בכלל המיקומים
ZH	须在医生处方下使用，不适用于所有地区
AR	يلزم وجود وصفة طبية، لا ينطبق في بعض المناطق الجغرافية
PT	Requer prescrição, não se aplica a todas as áreas geográficas
RU	Отпускается по рецепту (относится не ко всем регионам)

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Indications for use

The AutoSoft™ 30 infusion set is indicated for the subcutaneous infusion of insulin administered by an external pump. The infusion set is indicated for single use.

INTENDED USE

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

PRODUCT DESCRIPTION

AutoSoft 30 is an angled infusion set with an integrated insertion device. The built-in legs support a 30° insertion angle. This product is delivered ready for use. Both the insertion device and the separate tubing are sterile.

CONTRAINDICATIONS

AutoSoft 30 is neither intended nor indicated for use with blood, blood products or intravenous infusion (IV).

WARNINGS

- Carefully examine the instructions for use before inserting AutoSoft 30. Failure to follow instructions may result in pain or injury.
- AutoSoft 30 infusion sets with t:lock™ tubing connectors should only be used with Tandem Diabetes Care™ cartridges featuring the t:lock connector.
- AutoSoft 30 is only sterile and non-pyrogenic if the packaging is undamaged and unopened. Do not use it if the sterile paper or red seal has been broken.
- When using AutoSoft 30 for the first time, do so in the presence of a healthcare provider.
- Inaccurate medication delivery, site irritation or infection may result from improper insertion or maintenance of your infusion site.
- Do not re-insert the introducer needle into the infusion set. This could cause a tearing or damage of the soft cannula and may lead to unpredictable medication flow.
- Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.

- If the soft cannula bends during insertion, apply a new AutoSoft 30 immediately at a new site.

PRECAUTIONS

- Never try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.
- When priming the infusion set, make sure there are no air bubbles in the tubing. Consult your pump manufacturer's instructions for use.
- AutoSoft 30 is a single-use device and should be disposed of immediately after use. Do not clean or re-sterilize.
- Check the infusion set frequently to ensure that the soft cannula remains firmly in place. Replace with a new infusion set at a new site, if the infusion set is not in place. Since the cannula is soft, it will not cause any pain if it slips out, and this may happen without notice. The soft cannula must always be completely inserted to receive the full amount of medication.
- Do not put perfumes, deodorants, cosmetics or other products containing

disinfectants in contact with the connector or the tubing. This may impact the integrity of the infusion set.

- Replace the infusion set if the adhesive tape becomes loose or is dislodged from its original site.
- Replace the infusion set every two to three days, or per your healthcare provider's instructions.
- If your insertion site becomes irritated or infected, replace the infusion set and use a new site until the first site has healed.
- Wash your hands before disconnecting and reconnecting the infusion set. Consult your healthcare provider on how to compensate for missed medication while disconnected.
- Never point a loaded insertion device towards any body part, where insertion is not desired.

RECOMMENDATIONS

- Check your blood glucose level one to three hours after insertion of the AutoSoft 30. Measure your blood glucose on a regular basis. Discuss this with your healthcare provider.

- Do not change your infusion set just before bedtime, unless your blood glucose can be checked one to three hours after insertion.
- If your blood glucose level becomes unexplainably high or an occlusion alarm occurs, check for clogs and leaks. If in doubt, change your infusion set since the soft cannula could be dislodged. Discuss plan to bring down your blood glucose level with your healthcare provider should any of these problems arise. Test your blood glucose level to ensure problem is corrected.
- Carefully monitor your blood glucose level when disconnected and after reconnecting.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

STORAGE & DISPOSAL

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.

- Put the lid back on the insertion device and dispose in a proper sharp container according to local laws, to avoid risk of needle stick.
- Dispose the used infusion set according to local regulations for biohazardous waste.

Preparation and insertion

Read all instructions carefully before using AutoSoft 30. Consult your pump manufacturer's instructions for use before connecting AutoSoft 30 to your pump. Follow proper hygiene procedures. Please see the corresponding illustrations on pages 84-85.

CONTENTS

- A. Lid
- B. Needle protector
- C. Soft cannula
- D. Adhesive tape
- E. Cannula housing
- F. Connector
- G. Tubing
- H. Introducer needle
- I. Insertion device
- J. Sterile paper
- K. Lined indentations – hold here when loading AutoSoft 30 (see figure 7)
- L. Release buttons – press here when releasing AutoSoft 30 (see figure 9)
- M. Disconnect cover
- N. Circular protective cap

INSTRUCTIONS

Unfold the cover of this booklet to view images that correspond to the steps below.

1. Wash your hands before inserting AutoSoft 30.
2. Suggested insertion areas. Do not use an area right next to your previous site of insertion. Please consult your healthcare provider for recommended insertion site and rotation.
3. Clean the area of insertion with a disinfectant as directed by your healthcare provider. Be sure the area is air dry before inserting AutoSoft 30.
CAUTION: Check the tamper-proof seal and sterile paper. Do not use if the sterile paper or red seal has been broken.
4. Pull the red tape to remove the seal.
5. Pull off the sterile paper. Do not touch the spring at this stage.
6. Hold on the lined indentations with one hand and pull the lid off with the other.
CAUTION: Be careful not to bend or touch the AutoSoft 30 introducer needle.
IMPORTANT NOTE: If the infusion set is not

securely placed in the insertion device with the needle pointing straight ahead, pain or minor injury may occur during insertion.

7. Place your fingers on the lined indentations on each side. Pull back the spring softly until you hear a “click”. The needle protector will lift up or peel off. Do not touch or bend the introducer needle.
8. Immediately before insertion, carefully remove the needle protector. Be sure the tape is not stuck on the introducer needle. Check that the soft cannula does not exceed the introducer needle.
9. **Insertion:** Place your index finger on the top release button and your thumb on the bottom release button. Position AutoSoft 30 with its legs level (flat) on the skin to ensure an insertion angle of 30°. Make sure to maintain a 30° angle while inserting and gently squeeze once to insert AutoSoft 30.
10. Make sure that the cannula stays in place: Put your finger gently on the clear window while removing the insertion device and introducer needle. Pull it gently straight back.

11. Secure the cannula with one finger on the clear window. Carefully remove the backing paper under the cannula housing as shown.
12. Remove the backing paper under the front end. Gently massage the adhesive tape onto the skin.

IMPORTANT: Do not attempt to straighten the adhesive if curled. Start over with a new infusion set.

13. Put the lid back on for disposal.

14. **TUBING:** Remove paper from the tubing and the circular protective cap from the connector. Prime the tubing until medication dispenses from the tip. Ensure there are no air bubbles in the tubing. Consult your pump manufacturer's instructions for use for proper priming procedures.

15. Place a finger in front of the cannula housing and insert the connector needle. Make sure the connector locks with a "click". Fill the empty space in the cannula with the following prime:

U-100 INSULIN:
0.7 units (0.007 ml)

DISCONNECT

- A. AutoSoft 30 allows you to temporarily disconnect from your pump. Place a finger just in front of the cannula housing and gently squeeze the sides of the connector. Pull the connector needle straight out.
- B. Insert the circular protective cap into the connector needle until you hear a "click".
- C. Insert the disconnect cover into the cannula housing until you hear a "click".

RECONNECT

- D. Remove the circular protective cap. If necessary, prime the tubing until medication dispenses from the needle.
CAUTION: When priming, hold the tubing with the needle pointing down.
- E. Remove the disconnect cover from the cannula housing.
- F. Place a finger in front of the cannula housing. Push the connector needle straight in until you hear a "click".

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© 2023 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, AutoSoft, and t:lock are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. / Tous droits réservés. Tandem Diabetes Care, AutoSoft et t:lock sont des marques déposées ou des marques de commerce de Tandem Diabetes Care, Inc. aux États-Unis et/ou dans d'autres pays. AW-1000220_E
Tandem Diabetes Care, Inc.