



Ostomy

Digital Product Guide

Creating Healthy Bonds with Complete
Ostomy Products, Care and Support

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)



Helping to create healthy bonds

Together, we can help your patients create a healthy bond with their ostomy.

Convatec's skin-friendly and secure pouch connections empower your patients with the security to go out, live life, and connect with others. In addition to providing a comprehensive portfolio of ostomy product solutions, we are committed to elevating the experience of living with an ostomy through our me+™ program, which provides your patients with the education, support and resources they need.

Because, when your patient has a healthy bond with their ostomy, they have a healthy bond to their life.

LEARN MORE

to see how Convatec is creating Healthy Bonds



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

TABLE OF CONTENTS

me+™ SUPPORT 1

- me+™ Program Overview 3
- me+™ Virtual Support 5
- me+™ Style 7

EDUCATION 9

Ostomy Overview 11

- Introduction to Ostomy Care 11
- Choosing the Right Ostomy System 13

Guide to Ostomy Systems 15

- A Guide to Adhesives 15
- A Guide to Skin Barriers 19
- Moldable Skin Barriers 21
- Accordion Skin Barriers 23
- The Purpose of Convexity 25
- The Characteristics of Convexity 27
- Ostomy Pouch Options 29
- Common Pouch Features 31

TWO-PIECE SYSTEMS 33

- Natura™, Natura™+,
and Sur-FIT® Natura™ Systems 35
- Esteem Synergy® Systems 49

ONE-PIECE SYSTEMS 53

- Esteem Body™ Pouches 55
- Esteem™+ Pouches 61
- ActiveLife® Pouches 71

ACUTE KITS 75

PEDIATRIC PRODUCTS 81

ACCESSORIES 89

WOUND & FISTULA POUCHES 103

RESOURCES 107

- Patient Ostomy Resources & Support 109
- QR Code Quick Links 110
- Medicare Utilization Chart 111
- Product Code Index 113
- References 119

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

There is a Need for Ongoing Patient Support

Most common ostomy self-care challenges are related to leakage, skin issues, frequency of changes and difficulty caring for an ostomy¹

Following an interaction with a **me+™ WOC nurse:**

99% of survey respondents said they would recommend this service²

96% of those surveyed had a better understanding of how to care for their ostomy after a me+™ virtual visit (telehealth)³

87% of those surveyed said they avoided an additional visit to a healthcare facility²



THE me+™ PROGRAM SUPPORTS PATIENTS

The me+™ program is a comprehensive support program designed to support clinical care by providing patients with support, inspiration, insights and ideas throughout their entire ostomy journey.


me+™ program members can:

- ✓ Request free ostomy samples
- ✓ Speak to a WOC nurse or product specialist
- ✓ Receive complimentary virtual telehealth support
- ✓ Access online resources, education and guided recovery tools
- ✓ Connect to our Language Line Services, which offers translation support in over 240 foreign languages and dialects - at bedside or home

Health Care Professionals

TO ENROLL YOUR PATIENTS in the me+™ program

call us at 877-585-0470 or email concierge@Convatec.com



LEARN HOW

our me+™ program can provide the right ostomy support for your patient

TABLE OF CONTENTS

me+™ SUPPORT

me+™ SUPPORT

- me+™ Program Overview
- me+™ Virtual Support
- me+™ Style

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



Our goal is to **partner with you and your patients to deliver best-in-class care** throughout the patient journey:

✓ **Pre-Surgery / Surgery**

- Access online resources at mepluscare.com
- Request me+™ patient Starter Kit and education materials
- Enroll in me+ program
- Learn about the me+ recovery series, developed to aid rehabilitation after surgery

✓ **Weeks Immediately After Surgery**

- Consult with me+ WOC nurses and product specialists
- Personalized Starter Kit
- me+ welcome email and support
- me+ recovery phase 1 exercise booklet*

✓ **First Year & Beyond**

- Explore online resources and tools
- Get inspired by the me+ community
- Sample new products
- Continue me+ recovery series exercises*


Throughout their journey, your patients can call 1-800-422-8811 for live support from me+™ WOC nurses and product specialists

A me+™ Concierge team member can assist Healthcare Professionals with:

- ✓ Personalized support, including finding an ostomy supply retailer
- ✓ Enrolling your patients into the me+™ program
- ✓ Providing a me+ Starter Kit and free product samples for your patients

TO LEARN MORE contact your me+ concierge directly at 877-585-0470

Ostomy care has met its app



KEY FEATURES

- ✓ Nurse Login
- ✓ Product Selector
- ✓ Ship Samples
- ✓ Saved Kits
- ✓ me+™ Virtual Support
- ✓ Enroll Patients in me+™

Download the
Convatec Ostomy Nurse Solutions™ App Today



[OstomyNurseSolutions.com](https://www.OstomyNurseSolutions.com)

TABLE OF CONTENTS

me+™ SUPPORT

me+™ SUPPORT

- me+™ Program Overview
- me+™ Virtual Support
- me+™ Style

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

*The information in the series is not intended to provide medical advice. Patients and users should speak with their health care professional before engaging in physical activity following ostomy surgery. In the event they experience pain or discomfort, they should stop the activity immediately and contact their health care professional.



me+™ offers telehealth services nationwide

The me+™ program is the first of its kind to offer Telehealth in all 50 states and Puerto Rico. Our complimentary telehealth services provide visual support and live demonstrations through our telehealth partner, Corstrata™.

Interested in learning more about the me+™ virtual support option?



WATCH THE VIDEO

or call us today at 1-800-422-8811 to see if a telehealth consultation is right for your patient.


"It's truly allowed me to have a more positive outlook, which was certainly suffering before the visit."

- Michael, me+™ member and telehealth user

me+™ virtual support is not intended to be a substitute for professional medical advice, diagnosis, or treatment. Always seek the advice of your physician or other qualified health provider with any questions you may have regarding a medical condition.



Telehealth visits will be managed through Convatec partner, Corstrata™. Only visits scheduled through a Convatec referral are complimentary.



You Can Feel Confident Recommending me+™ Virtual Support Because It Is:

- ✓ Provided by Board-Certified Ostomy Nurses
- ✓ Available in all 50 states and Puerto Rico
- ✓ HIPAA Compliant, Private & Secure
- ✓ Available in English and Spanish

TABLE OF CONTENTS

me+™ SUPPORT

me+™ SUPPORT

- me+™ Program Overview
- me+™ Virtual Support
- me+™ Style

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



Look great and feel confident

Ostomysecrets® offers a full line of stylish and functional undergarments, swimwear and skin care designed to support your ostomy system.

Who says living with an ostomy can't be stylish, sleek and even sexy? With me+™ style, you don't have to change what you wear or hide your ostomy pouch under bulky, unflattering clothes. Ostomysecrets.com is a one-of-a-kind clothing site created by people with ostomies, for people with ostomies. From swimwear to women's thong underwear to men's athletic boxers, it's a full line of apparel you won't find anywhere else. Ostomysecrets.com gives you everything you need to conceal and support your ostomy pouch, without sacrificing your personal style.

HAVE QUESTIONS?

Call an Ostomysecrets® specialist at 1-877-613-6246, chat online at ostomysecrets.com or send us an email at info@ostomysecrets.com

SHOP NOW

and **SAVE 15% OFF YOUR FIRST ORDER** with code **STYLE15**.

TABLE OF CONTENTS

me+™ SUPPORT

me+™ SUPPORT

- me+™ Program Overview
- me+™ Virtual Support
- **me+™ Style**

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

EDUCATION

Feel more confident supporting your ostomy patients.

Take advantage of our educational resources to learn about ostomy basics, importance of skin care, and more to help you elevate the patient experience.

- ✓ An Introduction to Ostomy Care
- ✓ Choosing the Right Ostomy System
- ✓ A Guide to Adhesives
- ✓ A Guide to Skin Barriers
- ✓ Moldable Skin Barriers
- ✓ Accordion Skin Barriers
- ✓ The Purpose of Convexity
- ✓ The Characteristics of Convexity
- ✓ Ostomy Pouch Options
- ✓ Common Pouch Features

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYS

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRO

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

An Introduction to OSTOMY CARE

There are 3 primary types of ostomy surgeries. Understanding the fundamental difference is important to select the right product and provide the right support.



COLOSTOMY (large intestine)

A colostomy is a surgical procedure that creates an opening on the abdomen through which the large intestine is redirected to bypass diseased or damaged parts of the large intestine. The large intestine is brought to the surface of the skin, and sutured to the abdomen to form a new exit for the stool from the bowel.

Things to know about a Colostomy

- ✓ Primarily formed due to colon or rectal cancer or trauma
- ✓ Can be temporary or permanent
- ✓ Stool tends to be semi-formed to formed and predictable 1-2 times per day

ILEOSTOMY (small intestine)

An ileostomy is a surgical procedure that creates an opening on the abdomen through which the end of your small intestine is redirected to bypass diseased or damaged parts of the digestive system. The small intestine is brought to the surface of the skin, and sutured to your abdomen to form a new exit for waste from the bowels.

Things to know about an Ileostomy

- ✓ Primarily formed due to inflammatory bowel disease (Crohn's, Ulcerative Colitis)
- ✓ Can be temporary or permanent
- ✓ Stool is liquid or pasty and constant throughout the day
- ✓ The small intestine contains digestive enzymes and acids, and the output may cause skin irritation if it comes in contact with the skin

UROSTOMY (urinary tract)

A urostomy is a surgical procedure that creates an opening on the abdomen through which urine is redirected away from a bladder that is damaged or diseased. The most common type of urostomy is an ileal conduit. During this procedure the bladder is removed and a piece of small intestine (ileum) is separated from the rest of the bowel. This piece is used as a conduit for urine to pass out of the body. The ureters are attached to this new passageway and one end is brought up to an opening made on the abdomen.

Things to know about a Urostomy

- ✓ Primarily formed due to cancer or inflammation of the bladder
- ✓ Is almost always permanent
- ✓ Urine flows through the stoma as it is produced which means output is constant
- ✓ Urine may cause skin irritation if it comes in contact with the skin

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care

- Choosing the Right Ostomy System

TWO-PIECE SYSTEMS

- A Guide to Adhesives

- A Guide to Skin Barriers

ONE-PIECE SYSTEMS

- Moldable Skin Barriers

- Accordion Skin Barriers

- The Purpose of Convexity

ACUTE KITS

- The Characteristics of Convexity

- Ostomy Pouch Options

PEDIATRIC PRODUCTS

- Common Pouch Features

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Choosing the Right OSTOMY SYSTEM

The ostomy pouch collects stoma output. The pouch comes as part of a two-piece or one-piece system.



TWO-PIECE SYSTEM



The two-piece system includes a skin barrier and a separate pouch that securely join together. This option allows the pouch to be removed to empty or change, without removing the skin barrier as frequently.

Two-piece systems are designed for those

- ✓ who want to change their pouch more frequently than their skin barrier
- ✓ with an ileostomy who need to empty/change the pouch frequently
- ✓ who exercise, travel and play sports

ONE-PIECE SYSTEM



The skin barrier and pouch are provided as an integrated solution that cannot be separated. It offers one simple and convenient step to attach or remove, providing a fresh start every time. This option is thin and flexible, and designed for more frequent system changes.

One-piece systems are designed for those

- ✓ who require no added pressure on the abdomen caused when connecting a skin barrier and a pouch, especially post-operatively
- ✓ who prefer to change their ostomy system frequently
- ✓ who prefer a flexible and low profile ostomy system

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- **Choosing the Right Ostomy System**
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

A Guide to ADHESIVES

Trust the brand that's the gold standard in hydrocolloids. It's the same brand that invented them. Convatec revolutionized ostomy care with the development of Advanced Hydrocolloids. 40 years later, no one has more experience creating healthy bonds.



STOMAHESIVE® (standard wear)

Designed for people with mostly solid output, our unique Stomahesive® technology is a hydrocolloid adhesive that gently adheres to the skin, protecting it from the damaging effects of stoma effluent. We also offer a Modified Stomahesive® option, which features a touch of mineral oil that makes the barrier more flexible during wear, and removal easier on the skin.

They're designed to

- ✓ Offer reliable peristomal skin protection
- ✓ Adhere to dry and wet skin
- ✓ Be comfortable on the skin during wear and gentle to skin upon removal

DURAHESIVE® (extended wear)

Designed for people with mostly liquid output, our unique Durahesive® technology is a hydrocolloid adhesive that swells up upon hydration to protect the stoma. This special effect called turtlenecking, helps create a more secure seal without harming the skin. We also offer Durahesive® Plus, which features an enhanced tack, strength and cohesion designed for extended wear time or aggressive output.

They're designed to

- ✓ Be highly durable and resistant to liquid output
- ✓ Have increased adhesive strength
- ✓ Gently swell when in contact with liquid stool or urine to help create a more secure seal without harming the skin

TAPE COLLARS
are available in two forms:

Acrylic: Tape collar designed with good tack
Hydrocolloid: Designed for those with sensitivity to other tape collars

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- **A Guide to Adhesives**
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

A Guide to ADHESIVES

ABOUT OUR ADHESIVES

Our formulations are based on:

- ✓ Over four decades of proven and trusted ostomy solutions
- ✓ Leading adhesive technology, formulated after 15 years of research and with over 400 formulations
- ✓ Unmatched support designed to dramatically improve life with an ostomy



Choose
Durahesive®
if your output is
mostly liquid
stool or urine



Choose
Stomahesive®
if your output is
solid stool

[LEARN MORE](#)

about Convatec Adhesives



Adhesives Selection Guide

Convatec's proprietary adhesives are integral in all of our ostomy skin barriers. Our hydrocolloid adhesives are designed for use in high-moisture environments and have decades of proven product performance.

The chart below will help determine which adhesive may be best for your patient.

| Adhesive Type | Formed Stool | Semi-Formed Stool | Semi-Liquid Stool | Liquid Stool | Urine | System Type | |
|-----------------------|--------------|-------------------|-------------------|--------------|-------|---------------------|---------------|
| Durahesive® Plus | | | ✓ | ✓ | | 1-Piece Only | Extended Wear |
| Durahesive® | | ✓ | ✓ | ✓ | ✓ | 1-Piece and 2 Piece | |
| Stomahesive® | ✓ | ✓ | ✓ | ✓ | | 1-Piece and 2-Piece | Standard Wear |
| Modified Stomahesive® | ✓ | ✓ | | | | 1-Piece Only | |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- **A Guide to Adhesives**
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

A Guide to SKIN BARRIERS

The skin barrier is the piece of the pouching system that sticks to the body. It should fit snugly around the stoma to help prevent leakage and have an opening customized to stoma size and shape.



MOLDABLE



Proven to help patients break the cycle of leakage and skin damage.^{4,5,6} Convatec Moldable Technology™ is the only personalized skin barrier that offers a flexible adhesive that hugs the stoma contour, so it's easier to get a perfect, secure fit every time.

Designed for those

- ✓ who have large, round, oval or irregular shaped stomas
- ✓ with dexterity issues who find cutting challenging

CUT-TO-FIT




Cut-to-Fit barriers can be cut to match the shape and size of individual stomas. Using scissors, one can cut a circle, oval or any shape needed to fit the stoma.

Designed for those

- ✓ with a stoma size in-between pre-cut sizes
- ✓ with an irregular stoma shape

PRE-CUT



Pre-cut barriers have a uniform, circular stoma hole and are ready to apply right out of the box. They are available for the most common stoma sizes.

Designed for those

- ✓ with an evenly round stoma that fits the available pre-cut size
- ✓ with dexterity issues who find cutting challenging

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- **A Guide to Skin Barriers**
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYS

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRO

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

MOLDABLE

Skin Barriers

Convatec Moldable Technology™ is the only personalized skin barrier that offers a flexible adhesive that hugs the stoma contour, so it's easier to get a perfect, secure fit every time.

95% of new users rated it excellent for comfort, ease of preparing, attaching, removing, and reliability⁶

96% of people rated moldable barriers as "good" or "excellent" for ease of use⁶

CLINICALLY PROVEN

to help break the cycle of leakage and skin damage^{4,5,6}

As many as **80%** of individuals with an ostomy will experience peristomal skin complications⁷

- **95% of people** who started on Convatec Moldable Technology™ kept their skin healthy⁸
- **86% of people** who switched to Convatec Moldable Technology™ improved their skin condition⁸

Our Moldable products are highlighted with this colored bar throughout the Guide



Convatec's Unique Moldable Technology is made up of Tri-Laminate Protection



LEAK-LOCK™

First-of-its kind system provides a snug "turtleneck" effect for maximum leakage control and added security



PATENTED REBOUNDED MEMORY TECHNOLOGY™

Customizable so that it molds to the stoma for a perfect fit every time



STAYS ON

Clinically proven to prevent and improve peristomal skin irritation

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- **Moldable Skin Barriers**
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

See our Moldable Technology in action and

[REQUEST A FREE SAMPLE NOW](#)

[Request Product Sample](#)

ACCORDION

Skin Barriers

The Accordion Flange is an option on a range of Natura® skin barriers. The accordion flange lifts to provide generous finger room so coupling is easy and comfortable.



PRODUCT FEATURES

- 2x the lifting height of another lifting flange on the market?
- Simple to measure and cut or mold
- Range of flat, convex, Cut-to-Fit and Moldable options
- Compatible with SUR-FIT® Natura™ and Natura™+ pouch option

Designed for those

- ✓ with dexterity issues as it provides finger room so coupling is easy and comfortable
- ✓ to prevent pressure or trauma if abdominal tenderness is present due to post-operative stoma surgery/incisions

HOW IT WORKS

- Designed to help minimize abdominal pressure and pain when attaching the pouch.
- Connects with an audible 'click' and physical 'snap' to confirm the pouch is firmly attached.



Lift 'UP' flange to easily attach the pouch away from the body



Low profile in the 'DOWN' position so it's not visible under clothing

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- **Accordion Skin Barriers**
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYS

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRO

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

The Purpose of CONVEXITY



WHAT IS CONVEXITY?

All bodies are different. Convex wafers are designed to help manage irregular stomas and abdomens to help reduce leakage and simplify stoma management.

- Convexity is the outward curving of the baseplate or skin barrier.^{10,11}
- Convexity applies a uniform tension to the immediate peristomal skin, improving the projection of the stoma for drainage.¹²
- Convexity helps maintain a secure seal between the pouch and the skin.¹³

GROWING CONVEXITY NEED

Data shows the need for convex skin barriers growing year on year.¹⁴ This is due to a number of factors including:


- Changing surgical techniques/outcomes of stoma shape & profile
- Increased obesity rates

WHEN TO USE CONVEXITY



Flush Stomas

Located at the same level as the abdomen skin surface, and does not protrude¹⁵



Retracted Stomas

Stoma pulled into the abdomen below skin level or having little or no spout¹⁶



Skin Folds

A stoma which is deep in a skin fold/crease¹³



Soft or Flaccid Abdomens

An abdomen requiring firm level of support to stabilize the pouching system¹⁵

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- **The Purpose of Convexity**
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYS

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRO

ACCESSORIES


WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

The Characteristics of CONVEXITY

CHARACTERISTICS OF CONVEXITY¹⁷



Depth

The depth of the convex skin barrier is defined as the measurement from the apex of the dome to the base.




Compressibility

Compressibility is the capacity of the convex dome to be displaced or flattened.




Flexibility

Flexibility is how easily the convex skin barrier can bend.



Tension Location


Tension location is the position in which the convex dome exerts downward and outward forces on the peristomal topography.




Slope

The slope is the angle from the base of the convex skin barrier to the apex of the dome.

Tension location depths allow ostomy skin barriers to:



Fit a flush or retracted stoma to provide tension as close as possible to the stoma if protrusion of the stoma is needed.



Flatten peristomal skin folds using tension away from the stoma that may interfere with achieving a consistent, reliable seal.

Choosing the right convexity solution for different body profiles can be challenging. Please feel free to consult with one of our me+™ Specialists or WOC Nurses for further guidance.

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- Common Pouch Features

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Ostomy Pouch OPTIONS

Convatec offers a wide range of pouching options across both 2 Piece and 1 Piece Systems.



CLOSED



A closed pouch doesn't require draining. Once they are full, you simply remove and replace with a new pouch.

Additionally, ideal during physical activity, intimate moments, swimming, and travel.

Designed for those

- ✓ who have thick, formed stool that is predictable (1-2 times per day)

DRAINABLE



Drainable pouches allow you to drain the output periodically throughout the day so that you do not have to change pouches as frequently.

Designed for those

- ✓ who have loose, liquid stool or semi-formed stool that can be unpredictable or constant throughout the day

UROSTOMY



A Urostomy pouch has a faucet tap, soft tap, or fold over tap at the bottom that allows for draining urine as frequently as needed.

Designed for those

- ✓ with a Urostomy

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- **Ostomy Pouch Options**
- Common Pouch Features

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

Common

POUCH FEATURES

POUCH COLOR

Transparent (Clear front of pouch)

Clear pouch front for easy output assessment and determining when a pouch needs to be changed or emptied. Soft and quiet material on the body side of the pouch.

Opaque (Tan or warm gray)

Soft and quiet material cover both the front and the body side of the pouch for more discretion and comfort.



Easy-View Window

Some options come with an easy-view window for accurate positioning and observation.

FILTER OR NO FILTER

Filters help release and deodorize gas and keep the closed or drainable pouch flat. Note that the filter must be covered with a filter cover while swimming, bathing or showering.



DRAINABLE POUCH CLOSURE



InvisiClose™ Clipless Tail Closure

Self closing drainable pouch closure which provides an easy closure option that does not require a separate clip for added security.



Lock-it Pocket™

Folds inside out to hide away the secured drainage tail—providing the flexibility of a drainable pouch, with the feel and comfort of a closed-end pouch.

Traditional drainable pouches are closed with a separate clip.

UROSTOMY POUCH DRAINAGE

The Accuseal™ Tap with valve is a faucet type tap designed for easy manipulation and emptying. Other drainage options include a Soft Tap or Fold Over Tap.

All Urostomy pouches include an anti-reflux valve which is designed to help prevent backflow.



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

EDUCATION

- An Introduction to Ostomy Care
- Choosing the Right Ostomy System
- A Guide to Adhesives
- A Guide to Skin Barriers
- Moldable Skin Barriers
- Accordion Skin Barriers
- The Purpose of Convexity
- The Characteristics of Convexity
- Ostomy Pouch Options
- **Common Pouch Features**

TWO-PIECE SYS

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRO

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

TWO-PIECE

Systems



NATURA™

The Natura™ brand (inclusive of Natura™, Sur-Fit Natura™, and Natura™+) is an extensive range of skin barriers and pouches that connect together with a plastic ring called a flange in a simple snap-off, snap-on design.



Two-piece systems include a **skin barrier and a separate pouch** that securely join together. There are two brands that make up Convatec's Two-Piece Ostomy Systems: Natura™+ and Esteem Synergy®+.



ESTEEM SYNERGY®

The Esteem Synergy®+ brand are skin barriers and pouches that connect together with Adhesive Coupling Technology™. This is a flat, flexible adhesive foam ring on the pouch that sticks to a flat 'landing zone' on the skin barrier.

Color Key & Sizing Guide

| COLOR MATCH | | | | | | |
|-------------|------------------|------------------|------------------|------------------|------------------|---------------|
| FLANGE SIZE | 1 1/4" (32mm) | 1 1/2" (38mm) | 1 3/4" (45mm) | 2 1/4" (57mm) | 2 3/4" (70mm) | 4" (100mm) |

Color Key & Sizing Guide

| COLOR MATCH | | | |
|---------------|------------------------|-------------------------|------------------------|
| COUPLING SIZE | Small Stoma Opening | Medium Stoma Opening | Large Stoma Opening |



When ordering a Two-Piece System, **make sure you select the same flange/coupling size for the skin barrier and pouch.** Use the colored dots in the table to help.

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- Natura™+ Pouches
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

[Request Product Sample](#)

Natura™ Durahesive® Skin Barriers

Compatible with Natura™ +
and SUR-FIT® Natura™ pouches



SUR-FIT® NATURA™ SKIN BARRIER

Cut-to-Fit

- Durahesive® (Extended Wear)
- Quantity = 10 (*Box of 5)

HCPCS A4409

(413157 and 401905 are HCPCS A4410)

| Flange Size | Cutting Range | Acrylic Tape Collar WHITE | Acrylic Tape Collar TAN | Without Tape Collar |
|-----------------|---------------|---------------------------|-------------------------|---------------------|
| ● 1 1/4" (32mm) | 1/2 - 3/4" | 413159 | 413164 | 413153 |
| ● 1 1/2" (38mm) | 1/2 - 1" | 413160 | 413165 | 413154 |
| ● 1 3/4" (45mm) | 1/2 - 1 1/4" | 413161 | 413166 | 413155 |
| ● 2 1/4" (57mm) | 1/2 - 1 3/4" | 413162 | 413167 | 413156 |
| ● 2 3/4" (70mm) | 1/2 - 2 1/4" | 413163 | 413168 | 413157 |
| ● 4" (100mm) | 1/2 - 3 1/2" | N/A | N/A | 401905* |



NATURA™ SKIN BARRIER

Cut-to-Fit | Accordion Flange

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4409

| Flange Size | Cutting Range | Acrylic Tape Collar WHITE |
|-----------------|---------------------------|---------------------------|
| ● 1 3/4" (45mm) | 1/2 - 13/16" (13mm-21mm) | 421454 |
| ● 2 1/4" (57mm) | 1/2 - 1 5/16" (13mm-33mm) | 421458 |
| ● 2 3/4" (70mm) | 1/2 - 1 3/4" (13mm-45mm) | 421462 |



NATURA™ SKIN BARRIER

Moldable | Accordion Flange

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4409

| Flange Size | Stoma Size | Hydrocolloid Tape Collar TAN |
|-----------------|--------------------------|------------------------------|
| ● 2 1/4" (57mm) | 1/2 - 7/8" (13-22mm) | 421039 |
| ● 2 1/4" (57mm) | 7/8 - 1 1/4" (22-33mm) | 421040 |
| ● 2 3/4" (70mm) | 1 1/4 - 1 3/4" (33-45mm) | 421041 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- Natura™ + Pouches
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



NATURA™ SKIN BARRIER

**Cut-to-Fit | Convex |
Accordion Flange**

- Durahesive® (Extended Wear)
- Convex Wafer Depth = 6.65mm
- Quantity = 10

HCPCS A4407

| Flange Size | Cutting Range | Acrylic Tape Collar TAN | Hydrocolloid Tape Collar TAN |
|-----------------|------------------------|-------------------------------|------------------------------------|
| ● 1 3/4" (45mm) | 1/2 - 7/8" (13-22mm) | 421638 | 421639 |
| ● 2 1/4" (57mm) | 1/2 - 1 3/8" (13-35mm) | 421640 | 421641 |
| ● 2 3/4" (70mm) | 1/2 - 1 7/8" (13-48mm) | 421642 | 421643 |



NATURA™ SKIN BARRIER

Moldable

- Durahesive® (Extended Wear)
- Quantity = 10

HCPCS A4409

| Flange Size | Stoma Size | Acrylic Tape Collar WHITE | Hydrocolloid Tape Collar WHITE |
|-----------------|---|---------------------------------|--------------------------------------|
| ● 1 3/4" (45mm) | 1/2 - 7/8" (13-22mm) | 413417 | 411800 |
| ● 1 3/4" (45mm) | 7/8 - 1 1/4" (22-33mm) | 413418 | 411802 |
| ● 2 1/4" (57mm) | 1 1/4 - 1 3/4" (33-45mm) | 413419 | 411804 |
| ● 2 3/4" (70mm) | 1 1/4 - 1 3/4" by 1 3/4 - 2 1/8" (33x45mm by 45x56mm) Oval shape | 413420 | 411806 |



NATURA™ SKIN BARRIER

Moldable | Convex

- Durahesive® (Extended Wear)
- Quantity = 10


HCPCS A4407

| Flange Size | Stoma Size | Acrylic Tape Collar WHITE | Convex Wafer Depth |
|-----------------|--------------------------|---------------------------------|--------------------------|
| ● 1 3/4" (45mm) | 1/2 - 7/8" (13-22mm) | 404592 | 6.35mm |
| ● 1 3/4" (45mm) | 7/8 - 1 1/4" (22-33mm) | 404593 | 7mm |
| ● 2 1/4" (57mm) | 1 1/4 - 1 3/4" (33-45mm) | 404594 | 7mm |

SUR-FIT® NATURA™ SKIN BARRIER

Pre-cut | Convex

- Durahesive® (Extended Wear)
- Convex Wafer Depth = 6.65mm
- Quantity = 10



HCPCS A4407

| Flange Size | Stoma Size | Acrylic Tape Collar WHITE |
|-----------------|----------------|------------------------------|
| ● 1 3/4" (45mm) | 1/2" (13 mm) | 413177 |
| ● 1 3/4" (45mm) | 5/8" (16 mm) | 413178 |
| ● 1 3/4" (45mm) | 3/4" (19 mm) | 413179 |
| ● 1 3/4" (45mm) | 7/8" (22 mm) | 413180 |
| ● 1 3/4" (45mm) | 1" (25 mm) | 413181 |
| ● 1 3/4" (45mm) | 1 1/8" (28 mm) | 413182 |
| ● 1 3/4" (45mm) | 1 1/4" (32 mm) | 413183 |
| ● 1 3/4" (45mm) | 1 3/8" (35 mm) | 413184 |
| ● 2 1/4" (57mm) | 1 1/2" (38 mm) | 413185 |
| ● 2 1/4" (57mm) | 1 5/8" (41 mm) | 413186 |
| ● 2 1/4" (57mm) | 1 3/4" (45 mm) | 413187 |
| ● 2 1/4" (57mm) | 2" (50 mm) | 413188 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- Natura™ + Pouches
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Natura™ Stomahesive® Skin Barriers

Compatible with Natura™ +
and SUR-FIT® Natura™ pouches



SUR-FIT® NATURA™ SKIN BARRIER


Cut-to-Fit

- Stomahesive® (Standard Wear)
- Quantity = 10 (*Box of 5)

HPCPS A4414

(401577 and 401906 are HPCPS A4415)

| Flange Size | Cutting Range | Acrylic Tape Collar WHITE | Acrylic Tape Collar TAN | Without Tape Collar |
|-----------------|---------------|---------------------------|-------------------------|---------------------|
| ● 1 1/4" (32mm) | 1/2 - 3/4" | 125257 | 125262 | 401573 |
| ● 1 1/2" (38mm) | 1/2 - 1" | 125258 | 125263 | 401574 |
| ● 1 3/4" (45mm) | 1/2 - 1 1/4" | 125259 | 125264 | 401575 |
| ● 2 1/4" (57mm) | 1/2 - 1 3/4" | 125260 | 125265 | 401576 |
| ● 2 3/4" (70mm) | 1/2 - 2 1/4" | 125261 | 125266 | 401577 |
| ● 4" (100mm) | 1/2 - 3 1/2" | N/A | N/A | 401906* |



SUR-FIT® NATURA™ SKIN BARRIER

Pre-cut

- Stomahesive® (Standard Wear)
- Quantity = 10

HPCPS A4414

| Flange Size | Stoma Size | Acrylic Tape Collar TAN |
|-----------------|----------------|-------------------------|
| ● 1 3/4" (45mm) | 1/2" (13 mm) | 125267 |
| ● 1 3/4" (45mm) | 5/8" (16 mm) | 125268 |
| ● 1 3/4" (45mm) | 3/4" (19 mm) | 125269 |
| ● 1 3/4" (45mm) | 7/8" (22 mm) | 125270 |
| ● 1 3/4" (45mm) | 1" (25 mm) | 125271 |
| ● 1 3/4" (45mm) | 1 1/8" (28 mm) | 125272 |
| ● 1 3/4" (45mm) | 1 1/4" (32 mm) | 125273 |
| ● 1 3/4" (45mm) | 1 3/8" (35 mm) | 125274 |
| ● 2 1/4" (57mm) | 1 1/2" (38 mm) | 125275 |
| ● 2 1/4" (57mm) | 1 5/8" (41 mm) | 125276 |
| ● 2 1/4" (57mm) | 1 3/4" (45 mm) | 125277 |



NATURA™ SKIN BARRIER

Cut-to-Fit | Accordion Flange

- Stomahesive® (Standard Wear)
- Quantity = 10

HPCPS A4414

| Flange Size | Cutting Range | Acrylic Tape Collar WHITE |
|-----------------|-------------------------|---------------------------|
| ● 1 3/4" (45mm) | 1/2 - 13/16" (13-21mm) | 421456 |
| ● 2 1/4" (57mm) | 1/2 - 1 5/16" (13-33mm) | 421460 |
| ● 2 3/4" (70mm) | 1/2 - 1 3/4" (13-45mm) | 421464 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- **Natura™ Stomahesive® Skin Barriers**
- Natura™ + Pouches
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



NATURA™ SKIN BARRIER

Moldable | Accordion Flange

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

| Flange Size | Stoma Size | Hydrocolloid Tape Collar TAN |
|-----------------|--------------------------|---------------------------------|
| ● 2 1/4" (57mm) | 1/2 - 7/8" (13-22mm) | 421033 |
| ● 2 1/4" (57mm) | 7/8 - 1 1/4" (22-33mm) | 421034 |
| ● 2 3/4" (70mm) | 1 1/4 - 1 3/4" (33-45mm) | 421035 |




NATURA™ SKIN BARRIER

Moldable

- Stomahesive® (Standard Wear)
- Quantity = 10

HCPCS A4414

| Flange Size | Stoma Size | Acrylic Tape Collar WHITE | Hydrocolloid Tape Collar WHITE |
|-----------------|---|---------------------------------|--------------------------------------|
| ● 1 3/4" (45mm) | 1/2 - 7/8" (13-22mm) | 413421 | 411801 |
| ● 1 3/4" (45mm) | 7/8 - 1 1/4" (22-33mm) | 413422 | 411803 |
| ● 2 1/4" (57mm) | 1 1/4 - 1 3/4" (33-45mm) | 413423 | 411805 |
| ● 2 3/4" (70mm) | 1 1/4 - 1 3/4" by 1 3/4 - 2 1/8" (33x45mm by 45x56mm) Oval shape | 413424 | 411807 |



WATCH NOW

to learn how to apply
a Natura™ Moldable
Skin Barrier

Natura™ + Pouches

Compatible with SUR-FIT® Natura™
and Natura™ skin barriers



NATURA™ + POUCH

Drainable | InvisiClose™ Tail Closure

- Pouch Length = 12" (*pouch length = 14")
- Quantity = 10

HCPCS A5063 | No Filter

| | Flange Size | Transparent | Opaque |
|---|---------------|-------------|--------|
| ● | 1 3/4" (45mm) | 416418 | 421738 |
| ● | 2 1/4" (57mm) | 416421 | 421739 |
| ● | 2 3/4" (70mm) | 416424 | 421740 |
| ● | 4" (100mm) | 416472* | N/A |

HCPCS A4425 | Filter (Filter covers included)

| | Flange Size | Transparent | Opaque w/ Easy-View Window | Opaque |
|---|---------------|-------------|-------------------------------|--------|
| ● | 1 1/2" (38mm) | 416414 | 421746 | 416415 |
| ● | 1 3/4" (45mm) | 416416 | 421747 | 416417 |
| ● | 2 1/4" (57mm) | 416419 | 421748 | 416420 |
| ● | 2 3/4" (70mm) | 416422 | 421749 | 416423 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- **Natura™ Stomahesive® Skin Barriers**
- **Natura™ + Pouches**
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

SUR-FIT® NATURA™ POUCH

Drainable | Clipped Closure

- Pouch Length = 10"
- 2 Tail Closures per box
- Two-sided comfort panels in bullets
- Quantity = See table below



HCPCS A5063 | No Filter

| | Flange Size | Opaque Qty = 10 | Opaque Qty = 20 |
|---|---------------|-----------------|-----------------|
| ● | 1 1/4" (32mm) | 401505 | 404030 |
| ● | 1 1/2" (38mm) | 401506 | 404031 |
| ● | 1 3/4" (45mm) | 401507 | 404032 |
| ● | 2 1/4" (57mm) | 401508 | 404033 |
| ● | 2 3/4" (70mm) | 401509 | 404034 |

SUR-FIT® NATURA™ POUCH

Drainable | Clipped Closure

- Pouch Length = 12"
(*pouch length = 14" ; **pouch length = 10")
- 2 Tail Closures per box
- Quantity = See table below



HCPCS A5063 | No Filter

| | Flange Size | Transparent Qty = 10 | Opaque Qty = 10 | Opaque Qty = 20 | Opaque w/ 1 comfort panel Qty = 10 | Opaque w/ 1 comfort panel Qty = 20 |
|---|---------------|----------------------|-----------------|-----------------|------------------------------------|------------------------------------|
| ● | 1 1/4" (32mm) | 401510 | 401500 | 404025 | 401932 | 404014 |
| ● | 1 1/2" (38mm) | 401511 | 401501 | 404026 | 401933 | 404015 |
| ● | 1 3/4" (45mm) | 401512 | 401502 | 404027 | 401934 | 404016 |
| ● | 2 1/4" (57mm) | 401513 | 401503 | 404028 | 401935 | 404017 |
| ● | 2 3/4" (70mm) | 401514 | 401504 | 404029 | 401936 | 404018 |
| ● | 4" (100mm) | 401515* | N/A | N/A | N/A | N/A |

Continued ➤

HCPCS A4425 | Filter (Filter covers included)

| | Flange Size | Transparent Qty = 10 | Opaque Qty = 20 | Opaque Qty = 20** |
|---|---------------|----------------------|-----------------|-------------------|
| ● | 1 1/4" (32mm) | N/A | 411489 | 411288** |
| ● | 1 1/2" (38mm) | 411264 | 411490 | 411289** |
| ● | 1 3/4" (45mm) | 411265 | 411491 | 411290** |
| ● | 2 1/4" (57mm) | 411266 | 411492 | 411291** |
| ● | 2 3/4" (70mm) | 411267 | 411493 | 411292** |

SUR-FIT® NATURA™ POUCH

Drainable | High-Output

- Anti-Reflux Valve
- Outlet with spout and cap for high output drainage
- Pouch Length = 14"
- 1 Tail Closure, 1 Night Drainage Adapter, and 5 Replacement Filters per box
- Quantity = 5



HCPCS A4413 | Filter

| | Flange Size | Transparent | Opaque |
|---|---------------|-------------|--------|
| ● | 1 3/4" (45mm) | 420695 | 401557 |
| ● | 2 1/4" (57mm) | 420696 | 401558 |
| ● | 2 3/4" (70mm) | 420697 | 401559 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- **Natura™ + Pouches**
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



NATURA™ + POUCH

Closed-End

- Pouch Length = 8"
- Quantity = 30

HCPCS A5054 | No Filter

| | Flange Size | Opaque w/ Easy-View Window |
|---|---------------|----------------------------|
| ● | 1 1/4" (32mm) | 421677 |
| ● | 1 1/2" (38mm) | 421679 |
| ● | 1 3/4" (45mm) | 421681 |
| ● | 2 1/4" (57mm) | 421683 |
| ● | 2 3/4" (70mm) | 421893 |

HCPCS A4419 | Filter (Filter covers included)

| | Flange Size | Opaque w/ Easy-View Window | Opaque |
|---|---------------|-------------------------------|--------|
| ● | 1 1/4" (32mm) | 421676 | N/A |
| ● | 1 1/2" (38mm) | 421678 | 416403 |
| ● | 1 3/4" (45mm) | 421680 | 416406 |
| ● | 2 1/4" (57mm) | 421682 | 416409 |
| ● | 2 3/4" (70mm) | 421798 | 416412 |

NATURA™ + POUCH

Closed-End

- Pouch Length = 6"
- Quantity = 30

HCPCS A4419 | Filter (Filter covers included)

| | Flange Size | Opaque w/ Easy-View Window |
|---|---------------|----------------------------|
| ● | 1 3/4" (45mm) | 421796 |
| ● | 2 1/4" (57mm) | 421797 |

SUR-FIT® NATURA™ POUCH

Closed-End | Flange Cap

- Quantity = 25

HCPCS A5055 | Filter (Filter covers included)

| | Flange Size | Opaque w/ 1 comfort panel |
|---|---------------|---------------------------|
| ● | 1 3/4" (45mm) | 401909 |



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- **Natura™ + Pouches**
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



SUR-FIT® NATURA™ POUCH

Urostomy | Accuseal™ Tap with Valve

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10

HPCPS A4432

| | Flange Size | Transparent | Opaque | Opaque* |
|---|---------------|-------------|--------|---------|
| ● | 1 1/4" (32mm) | 401542 | 401551 | 401547* |
| ● | 1 1/2" (38mm) | 401543 | 401552 | 401548* |
| ● | 1 3/4" (45mm) | 401544 | 401553 | 401549* |
| ● | 2 1/4" (57mm) | 401545 | 401554 | 401550* |
| ● | 2 3/4" (70mm) | 401546 | 401555 | N/A |

NATURA™+ POUCH

Urostomy | Soft Tap

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Soft Tap Drainage Adapters per box
- Quantity = 10



HPCPS A5073

| | Flange Size | Transparent | Transparent* |
|---|---------------|-------------|--------------|
| ● | 1 1/4" (32mm) | N/A | 421074* |
| ● | 1 3/4" (45mm) | 413437 | N/A |
| ● | 2 1/4" (57mm) | 413438 | N/A |
| ● | 2 3/4" (70mm) | 413439 | N/A |

SUR-FIT® NATURA™ POUCH

Urostomy | Fold Over Tab

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Universal Night Drainage Adapters per box
- Quantity = 10



HPCPS A5073

| | Flange Size | Transparent | Transparent* |
|---|---------------|-------------|--------------|
| ● | 1 1/4" (32mm) | 401533 | 401538* |
| ● | 1 1/2" (38mm) | 401534 | 401539* |
| ● | 1 3/4" (45mm) | 401535 | 401540* |
| ● | 2 1/4" (57mm) | 401536 | 401541* |
| ● | 2 3/4" (70mm) | 401537 | N/A |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- **Natura™ + Pouches**
- Esteem Synergy® Skin Barriers
- Esteem Synergy® Pouches

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES


RESOURCES

Request Product Sample

Esteem Synergy[®] Skin Barriers

Compatible with Esteem Synergy[®] + Pouches

ESTEEM SYNERGY[®] SKIN BARRIER



Cut-to-Fit

- Stomahesive[®] (Standard Wear)
- Quantity = 10

HPCPS A4414

| Coupling Size | Cutting Range | Acrylic Tape Collar WHITE |
|------------------------|---------------------|------------------------------|
| ● Small Stoma Opening | up to 1 3/8" (35mm) | 405456 |
| ● Medium Stoma Opening | up to 1 7/8" (48mm) | 405457 |
| ● Large Stoma Opening | up to 2 3/8" (61mm) | 405458 |

HPCPS A4415

| Coupling Size | Cutting Range | Without Tape Collar |
|------------------------|---------------------|---------------------|
| ● Small Stoma Opening | up to 1 3/8" (35mm) | 405467* |
| ● Medium Stoma Opening | up to 1 7/8" (48mm) | 405468 |
| ● Large Stoma Opening | up to 2 3/8" (61mm) | 403948 |

*HPCPS is A4414

ESTEEM SYNERGY[®] SKIN BARRIER



Moldable

- Stomahesive[®] (Standard Wear)
- Quantity = 10

HPCPS A4414

| Coupling Size | Stoma Size | Hydrocolloid Tape Collar WHITE |
|------------------------|--------------------------|--------------------------------------|
| ● Medium Stoma Opening | 1/2 - 7/8" (13-22mm) | 411657 |
| ● Medium Stoma Opening | 7/8 - 1 1/4" (22-33mm) | 411658 |
| ● Large Stoma Opening | 1 1/4 - 1 3/4" (33-45mm) | 411659 |

ESTEEM SYNERGY[®] SKIN BARRIER




Moldable | Convex

- Durahesive[®] (Extended Wear)
- Convex Wafer Depth = 6.65mm
- Quantity = 10

HPCPS A4407

| Coupling Size | Stoma Size | Acrylic Tape Collar WHITE |
|------------------------|--------------------------|------------------------------|
| ● Small Stoma Opening | 1/2 - 7/8" (13-22mm) | 409268 |
| ● Medium Stoma Opening | 7/8 - 1 1/4" (22-33mm) | 409269 |
| ● Large Stoma Opening | 1 1/4 - 1 3/4" (33-45mm) | 409270 |

ESTEEM SYNERGY[®] SKIN BARRIER



Pre-Cut

- Stomahesive[®] (Standard Wear)
- Quantity = 10

HPCPS A4414

| Coupling Size | Stoma Size | Acrylic Tape Collar WHITE |
|------------------------|----------------|------------------------------|
| ● Small Stoma Opening | 5/8" (16 mm) | 405473 |
| ● Small Stoma Opening | 3/4" (19 mm) | 405474 |
| ● Small Stoma Opening | 7/8" (22 mm) | 405475 |
| ● Small Stoma Opening | 1" (25 mm) | 405476 |
| ● Small Stoma Opening | 1 1/8" (28 mm) | 405477 |
| ● Small Stoma Opening | 1 1/4" (32 mm) | 405478 |
| ● Small Stoma Opening | 1 3/8" (35 mm) | 405479 |
| ● Medium Stoma Opening | 1 1/2" (38 mm) | 405480 |
| ● Medium Stoma Opening | 1 5/8" (41 mm) | 405481 |
| ● Medium Stoma Opening | 1 3/4" (45 mm) | 405482 |

TABLE OF CONTENTS

me+[™] SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura[™] Durahesive[®] Skin Barriers
- Natura[™] Stomahesive[®] Skin Barriers
- Natura[™] + Pouches
- Esteem Synergy[®] Skin Barriers
- Esteem Synergy[®] Pouches

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Esteem Synergy® Pouches

Compatible with Esteem Synergy® Skin Barriers

ESTEEM SYNERGY® + POUCH

Drainable | InvisiClose™ Tail Closure

- Pouch Length = 12"
- Quantity = 10



HCPCS A5063 | No Filter

| Coupling Size | Transparent |
|-----------------------|-------------|
| ● Large Stoma Opening | 416903 |

HCPCS A4425 | Filter (Filter covers included)

| Coupling Size | Transparent | Opaque |
|------------------------|-------------|--------|
| ● Small Stoma Opening | 416795 | 416796 |
| ● Medium Stoma Opening | 416798 | 416799 |
| ● Large Stoma Opening | 416901 | 416902 |

ESTEEM SYNERGY® + POUCH

Closed-End

- Pouch Length = 8"
- Quantity = 30



HCPCS A4419 | Filter (Filter covers included)

| Coupling Size | Transparent | Opaque |
|------------------------|-------------|--------|
| ● Small Stoma Opening | 416774 | 416775 |
| ● Medium Stoma Opening | 416778 | 416779 |
| ● Large Stoma Opening | 416782 | 416783 |

ESTEEM SYNERGY® + POUCH

Closed-End

- Pouch Length = 6"
- Quantity = 30



HCPCS A4419 | Filter (Filter covers included)

| Coupling Size | Transparent | Opaque |
|------------------------|-------------|--------|
| ● Small Stoma Opening | 416776 | 416777 |
| ● Medium Stoma Opening | N/A | 416781 |

ESTEEM SYNERGY® POUCH

Urostomy | Accuseal™ Tap with Valve

- Anti-Reflux Valve
- Pouch Length = 10" (*pouch length = 9")
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10



HCPCS A4432

| Coupling Size | Transparent | Transparent* | Opaque |
|------------------------|-------------|--------------|--------|
| ● Small Stoma Opening | 405451 | 405446* | 405448 |
| ● Medium Stoma Opening | 405452 | 405447* | 405449 |
| ● Large Stoma Opening | 405453 | N/A | 405450 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE
SYSTEMS

TWO-PIECE SYSTEMS

- Natura™ Durahesive® Skin Barriers
- Natura™ Stomahesive® Skin Barriers
- Natura™ + Pouches
- Esteem Synergy® Skin Barriers
- **Esteem Synergy® Pouches**

ONE-PIECE SYS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

ONE-PIECE

Systems



EsteemBody™

with Leak Defense™

Introducing new Esteem Body™ soft convex one-piece ostomy systems with Leak Defense™. Our gold-standard adhesives are combined with a comprehensive soft convexity range designed to adapt to the body for a secure, longer-lasting seal. Our new modern 8-shaped pouch was designed with patients & nurses to sag less and provide more discretion.



ESTEEM™ +

The Esteem™ + brand is an all-in-one solution combining the skin barrier and pouch in a single unit. Designed to be easy and convenient, it's one simple step to attach or remove. Available in multiple adhesive options, growing convexity options, and with Convatec's latest pouch technologies.



ActiveLife®

Active Life® is a legacy one-piece system brand that also combines the skin barrier and pouch in a single unit.

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ONE-PIECE SYSTEMS

- Esteem Body™ Pouches
- Esteem™+ Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

EsteemBody™ Pouches



ESTEEM BODY™ POUCH

Cut-to-Fit | Drainable | Convex | InvisiClose™ Tail Closure

- Durahesive® (Extended Wear)
- Discreet warm gray color
- Convex Wafer Depth = 3.5 mm
- Pouch Length = 11" (286 mm)
- Pouch Capacity = 570 ml
- Belt Tabs
- Without Tape Collar
- Quantity = 10

HPCPS A5057 | Filter (Filter covers included)


| Cutting Range | Tension Location | Transparent | Opaque w/ Easy-View Window |
|------------------------|------------------|-------------|----------------------------|
| 3/8 - 2 1/8" (10-55mm) | 2 3/8" (60mm) | 423651 | N/A |
| 3/8 - 1 3/4" (10-45mm) | 2" (50mm) | 423652 | 423643 |
| 3/8 - 1 3/8" (10-35mm) | 1 1/2" (40mm) | 423653 | 423644 |
| 3/8 - 1" (10-25mm) | 1 1/8" (30mm) | 423654 | 423645 |

Pre-cut | Drainable | Convex

HPCPS A5057 | Filter (Filter covers included)

| Stoma Size | Tension Location | Opaque w/ Easy-View Window |
|---------------|------------------|----------------------------|
| 1" (25mm) | 1 1/2" (40mm) | 423646 |
| 1 1/8" (30mm) | 1 1/2" (40mm) | 423647 |
| 1 3/8" (35mm) | 1 1/2" (40mm) | 423648 |

Sizing Guidance for Esteem Body™ Pouches



Tension Location



ESTEEM BODY™ POUCH

Cut-to-Fit | Drainable | Convex | InvisiClose™ Tail Closure

- Durahesive® (Extended Wear)
- Discreet warm gray color
- Convex Wafer Depth = 7 mm
- Pouch Length = 11" (286 mm)
- Pouch Capacity = 570 ml
- Belt Tabs
- Without Tape Collar
- Quantity = 10

HPCPS A5057 | Filter (Filter covers included)

| Cutting Range | Tension Location | Transparent | Opaque w/ Easy-View Window |
|------------------------|------------------|-------------|----------------------------|
| 3/8 - 2 1/8" (10-55mm) | 2 3/8" (60mm) | 423664 | N/A |
| 3/8 - 1 3/4" (10-45mm) | 2" (50mm) | 423665 | 423656 |
| 3/8 - 1 3/8" (10-35mm) | 1 1/2" (40mm) | 423666 | 423657 |
| 3/8 - 1" (10-25mm) | 1 1/8" (30mm) | 423667 | 423658 |

Pre-cut | Drainable | Convex

HPCPS A5057 | Filter (Filter covers included)

| Stoma Size | Tension Location | Opaque w/ Easy-View Window |
|---------------|------------------|----------------------------|
| 1" (25mm) | 1 1/2" (40mm) | 423659 |
| 1 1/8" (30mm) | 1 1/2" (40mm) | 423660 |
| 1 3/8" (35mm) | 1 1/2" (40mm) | 423661 |

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ONE-PIECE SYSTEMS

- Esteem Body™ Pouches
- Esteem™+ Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



ESTEEM BODY™ POUCH

Cut-to-Fit | Closed-End | Convex

- Modified Stomahesive® (Standard wear)
- Discreet warm gray color
- Convex Wafer Depth = 3.5 mm
- Pouch Length = 9" (224 mm)
- Pouch Capacity = 520 ml
- Belt Tabs
- Without Tape Collar
- Quantity = 30

HCPCS A4417 | Filter (Filter covers included)

| Cutting Range | Tension Location | Opaque w/ Easy-View Window |
|------------------------|------------------|-------------------------------|
| 3/8 - 2 1/8" (10-55mm) | 2 3/8" (60mm) | 423676 |
| 3/8 - 1 3/4" (10-45mm) | 2" (50mm) | 423677 |
| 3/8 - 1 3/8" (10-35mm) | 1 1/2" (40mm) | 423678 |



ESTEEM BODY™ POUCH

Cut-to-Fit | Closed-End | Convex

- Modified Stomahesive® (Standard wear)
- Discreet warm gray color
- Convex Wafer Depth = 7 mm
- Pouch Length = 9" (224 mm)
- Pouch Capacity = 520 ml
- Belt Tabs
- Without Tape Collar
- Quantity = 30

HCPCS A4417 | Filter (Filter covers included)

| Cutting Range | Tension Location | Opaque w/ Easy-View Window |
|------------------------|------------------|-------------------------------|
| 3/8 - 2 1/8" (10-55mm) | 2 3/8" (60mm) | 423689 |
| 3/8 - 1 3/4" (10-45mm) | 2" (50mm) | 423690 |
| 3/8 - 1 3/8" (10-35mm) | 1 1/2" (40mm) | 423691 |



ESTEEM BODY™ POUCH

Cut-to-Fit | Urostomy

- Durahesive® (Extended Wear)
- Discreet warm gray color
- Convex Wafer Depth = 3.5 mm
- Pouch Length = 10" (259 mm)
- Pouch Capacity = 460 ml
- Belt Tabs
- Soft Tap Closure
- Without Tape Collar
- Quantity = 10

HCPCS A4393

| Cutting Range | Tension Location | Transparent | Opaque w/ Easy-View Window |
|------------------------|------------------|-------------|-------------------------------|
| 3/8 - 1 3/4" (10-45mm) | 2" (50mm) | 423712 | N/A |
| 3/8 - 1 3/8" (10-35mm) | 1 1/2" (40mm) | 423713 | 423704 |
| 3/8 - 1" (10-25mm) | 1 1/8" (30mm) | 423714 | 423705 |

Pre-cut | Urostomy | Convex**HCPCS A4393**

| Stoma Size | Tension Location | Transparent |
|---------------|------------------|-------------|
| 5/8" (15mm) | 1 1/8" (30mm) | 423715 |
| 3/4" (20mm) | 1 1/8" (30mm) | 423716 |
| 1" (25mm) | 1 1/8" (30mm) | 423717 |
| 1 1/8" (30mm) | 1 1/2" (40mm) | 423718 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE
SYSTEMS**ONE-PIECE SYSTEMS**

- Esteem Body™ Pouches
- Esteem™+ Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



ESTEEM BODY™ POUCH

Cut-to-Fit | Urostomy

- Durahesive® (Extended Wear)
- Discreet warm gray color
- Convex Wafer Depth = 7 mm
- Pouch Length = 10" (259 mm)
- Pouch Capacity = 460 ml
- Belt Tabs
- Soft Tap Closure
- Without Tape Collar
- Quantity = 10

HCPs A4393

| Cutting Range | Tension Location | Transparent | Opaque w/ Easy-View Window |
|------------------------|------------------|-------------|----------------------------|
| 3/8 - 1 3/4" (10-45mm) | 2" (50mm) | 423725 | N/A |
| 3/8 - 1 3/8" (10-35mm) | 1 1/2" (40mm) | 423726 | 423722 |
| 3/8 - 1" (10-25mm) | 1 1/8" (30mm) | 423727 | 423723 |



Request a **FREE** sample of
EsteemBody™
 with Leak Defense™

LEARN MORE

visit convatec.com
 call 877-585-0470 or email us at
concierge@Convatec.com

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

**ONE-PIECE
SYSTEMS**

ONE-PIECE SYSTEMS

- Esteem Body™ Pouches
- Esteem™+ Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Esteem™ + Pouches



ESTEEM™ + POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Durahesive® PLUS (Extended Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- Includes Belt Tabs
- Quantity = 10

HPCPS A5056 | Filter (Filter covers included)

| Cutting Range | Transparent | Opaque |
|------------------------|-------------|--------|
| 3/4 - 2 1/2" (19-64mm) | 416976 | 416975 |

Pre-cut | Drainable | InvisiClose™ Tail Closure

HPCPS A4424 | Filter (Filter covers included)

| Stoma Size | Opaque |
|---------------|--------|
| 1" (25mm) | 421025 |
| 1 1/8" (28mm) | 421026 |
| 1 1/4" (32mm) | 421027 |
| 1 3/8" (35mm) | 421028 |
| 1 1/2" (38mm) | 421029 |



ESTEEM™ + POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 12"
- Quantity = 10

HPCPS A5061 | No Filter

| Cutting Range | Transparent | Opaque w/ Easy-View Window |
|--------------------------|-------------|----------------------------|
| 13/16 - 2 3/4" (20-70mm) | 416721 | 421853 |

HPCPS A4424 | Filter

| Cutting Range | Transparent | Opaque w/ Easy-View Window | Opaque |
|--------------------------|-------------|----------------------------|--------|
| 13/16 - 2 3/4" (20-70mm) | 416718 | 421858 | 416719 |



ESTEEM™ + POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)*
- Without Tape Collar
- Pouch Length = 14"
- Quantity = 10

HPCPS A5061 | No Filter

| Cutting Range | Transparent |
|------------------|-------------|
| Up to 4" (100mm) | 416908 |

*Double-layer adhesive (Stomahesive body-side, Durahesive top-side) which combines the skin friendly aspects of Stomahesive with the turtlenecking protection of Durahesive.



ESTEEM™ + POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 12"
- 2 Tail Closures per box
- Quantity = 10

HPCPS A4424 | Filter (Filter covers included)*

| Cutting Range | Transparent |
|--------------------------|-------------|
| 13/16 - 2 3/4" (20-70mm) | 421848 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ONE-PIECE SYSTEMS

- Esteem Body™ Pouches
- Esteem™ + Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



ESTEEM™+ POUCH

**Pre-Cut | Drainable |
InvisiClose™ Tail Closure**

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 12"
- Quantity = 10

HPCPS A5061 | No Filter

| Stoma Size | Opaque w/ Easy-View Window |
|----------------|----------------------------|
| 1" (25mm) | 421854 |
| 1 3/16" (30mm) | 421855 |
| 1 3/8" (35mm) | 421856 |
| 1 9/16" (40mm) | 421857 |

HPCPS A4424 | Filter (Filter covers included)

| Stoma Size | Transparent | Opaque w/ Easy-View Window | Opaque |
|----------------|-------------|-------------------------------|--------|
| 1" (25mm) | 416724 | 421859 | 416725 |
| 1 3/16" (30mm) | 416728 | 421860 | 416729 |
| 1 3/8" (35mm) | 416732 | 421861 | 416733 |
| 1 9/16" (40mm) | 416736 | 421862 | 416737 |



ESTEEM™+ POUCH

**Pre-Cut | Drainable | Convex |
InvisiClose™ Tail Closure**

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Convex Wafer Depth = 5mm
- Pouch Length = 12"
- Includes Belt Tabs
- Quantity = 10

HPCPS A5057 | Filter (Filter covers included)

| Stoma Size | Transparent | Opaque |
|---------------|-------------|--------|
| 3/4" (19mm) | 416738 | 416739 |
| 7/8" (22mm) | 416740 | 416741 |
| 1" (25mm) | 416742 | 416743 |
| 1 1/8" (28mm) | 416744 | 416745 |
| 1 1/4" (32mm) | 416746 | 416747 |
| 1 3/8" (35mm) | 416748 | 416749 |
| 1 1/2" (38mm) | 416750 | 416751 |



ESTEEM™+ POUCH

**Moldable | Drainable |
InvisiClose™ Tail Closure**

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 12"
- Quantity = 10

HPCPS A5061 | No Filter

| Stoma Size | Transparent |
|-----------------------------|-------------|
| 13/16 - 1 3/16" (20-30mm) | 413518 |
| 1 3/16 - 1 9/16" (30-40mm) | 413519 |
| 1 9/16 - 1 15/16" (40-50mm) | 413520 |

HPCPS A4424 | Filter (Filter covers included)

| Stoma Size | Transparent | Opaque w/ Easy-View Window |
|-----------------------------|-------------|-------------------------------|
| 13/16 - 1 3/16" (20-30mm) | 413515 | 413521 |
| 1 3/16 - 1 9/16" (30-40mm) | 413516 | 413522 |
| 1 9/16 - 1 15/16" (40-50mm) | 413517 | 413523 |

**WATCH NOW**

to learn how to apply
an Esteem™+ One-Piece
Moldable pouch

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

**ONE-PIECE
SYSTEMS****ONE-PIECE SYSTEMS**

- Esteem Body™ Pouches
- **Esteem™+ Pouches**
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



ESTEEM™ + FLEX CONVEX POUCH

Cut-to-Fit | Drainable | Convex

- Esteem™ + Flex Convex Adhesive (Standard Wear)
- Without Tape Collar
- Convex Wafer Depth = 7mm
- Safe Seal Clipless Closure (1 Tail Closure per box)
- Pouch Length = 8 2/3"
- Includes Belt Tabs
- Quantity = 10


HCPCS A4389 | Filter (Filter covers included)

| Cutting Range | Plateau Size | Opaque w/ Easy-View Window |
|-------------------------|--------------|----------------------------|
| 3/4 - 1" (20-25mm) | V3 | 421622 |
| 3/4 - 1 3/8" (20-35mm) | V2 | 421619 |
| 3/4 - 1 2/3" (20-43mm) | V1 | 421615 |
| 3/4 - 2 6/25" (20-57mm) | V0 | 422159 |

Pre-Cut | Drainable | Convex

HCPCS A4389 | Filter (Filter covers included)

| Stoma Size | Plateau Size | Opaque w/ Easy-View Window |
|---------------|--------------|----------------------------|
| 1" (25mm) | V2 | 421620 |
| 1 1/6" (30mm) | V2 | 421621 |
| 1" (25mm) | V1 | 421616 |
| 1 1/6" (30mm) | V1 | 421617 |
| 1 3/8" (35mm) | V1 | 421618 |

Sizing Guidance for
Esteem™ + Flex and Soft Convex

ESTEEM™ + SOFT CONVEX POUCH

Cut-to-Fit | Drainable | Convex

- Esteem™ + Soft Convex Adhesive (Standard Wear)
- Without Tape Collar
- Convex Wafer Depth = 4mm
- Safe Seal Clipless Closure (1 Tail Closure per box)
- Pouch Length = 8 2/3"
- Includes Belt Tabs
- Quantity = 10

HCPCS A4389 | Filter (Filter covers included)

| Cutting Range | Plateau Size | Opaque w/ Easy-View Window |
|------------------------|--------------|----------------------------|
| 3/8 - 1 1/8" (10-28mm) | V3 | 422365 |
| 3/5 - 1 3/5" (15-40mm) | V2 | 422364 |
| 3/4 - 1 7/8" (20-47mm) | V1 | 422363 |

Pre-Cut | Drainable | Convex

HCPCS A4389 | Filter (Filter covers included)

| Stoma Size | Plateau Size | Opaque w/ Easy-View Window |
|---------------|--------------|----------------------------|
| 1" (25mm) | V2 | 422545 |
| 1 1/6" (30mm) | V2 | 422546 |
| 1" (25mm) | V1 | 422539 |
| 1 1/6" (30mm) | V1 | 422540 |
| 1 3/8" (35mm) | V1 | 422541 |

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE
SYSTEMS**ONE-PIECE SYSTEMS**

- Esteem Body™ Pouches
- Esteem™ + Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample



ESTEEM™+ POUCH

Cut-to-Fit | Closed-End

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 8"
- Quantity = 30

HPCPS A4416 | Filter (Filter covers included)

| Cutting Range | Transparent | Opaque w/ Easy-View Window | Opaque |
|-----------------------------|-------------|-------------------------------|--------|
| 13/16 - 2 3/4" (20-70mm) | 416700 | 421821 | 416701 |

ESTEEM™+ POUCH

Cut-to-Fit | Closed-End

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 6"
- Quantity = 30

HPCPS A4416 | Filter (Filter covers included)

| Cutting Range | Opaque w/ Easy-View Window |
|--------------------------|----------------------------|
| 13/16 - 2 3/4" (20-70mm) | 421816 |



ESTEEM™+ POUCH

Pre-Cut | Closed-End

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 8"
- Quantity = 30

HPCPS A4416 | Filter (Filter covers included)

| Stoma Size | Transparent | Opaque w/ Easy-View Window | Opaque |
|----------------|-------------|-------------------------------|--------|
| 1" (25mm) | 416703 | 421687 | 416704 |
| 1 3/16" (30mm) | 416706 | 421688 | 416707 |
| 1 3/8" (35mm) | 416709 | 421689 | 416710 |
| 1 9/16" (40mm) | 416712 | 421822 | 416713 |
| 2" (50mm) | N/A | 421823 | 416715 |



ESTEEM™+ POUCH

Moldable | Closed-End

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 8"
- Quantity = 30

HPCPS A4416 | Filter (Filter covers included)

| Stoma Size | Opaque w/ Easy-View Window |
|----------------------------|----------------------------|
| 13/16 - 1 3/16" (20-30mm) | 413509 |
| 1 3/16 - 1 9/16" (30-40mm) | 413510 |



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

**ONE-PIECE
SYSTEMS****ONE-PIECE SYSTEMS**

- Esteem Body™ Pouches
- **Esteem™+ Pouches**
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

**ESTEEM™ + POUCH****Urostomy | Cut-to-Fit | Accuseal™ Tap with Valve**

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Anti-Reflux Valve
- Pouch Length = 10"
- Includes Belt Tabs
- 2 Accuseal Night Drainage Adapters per box
- Quantity = 10

HPCPS A4428**ESTEEM™ + FLEX CONVEX POUCH****Urostomy | Cut-to-Fit | Convex**

- Esteem™ + Flex Convex Adhesive (Standard Wear)
- Without Tape Collar
- Convex Wafer Depth = 7mm
- Soft Tap Closure
- Pouch Length = 7 2/3"
- Includes Belt Tabs
- 2 Esteem™ + Flex Convex Night Drainage Adapters per box
- Quantity = 10

HPCPS A4392

| Cutting Range | Plateau Size | Opaque w/ Easy-View Window |
|------------------------|--------------|----------------------------|
| 3/8 - 1" (10-25mm) | V3 | 421628 |
| 3/8 - 1 3/8" (10-35mm) | V2 | 421626 |
| 3/8 - 1 2/3" (10-43mm) | V1 | 421623 |

Urostomy | Pre-Cut | Convex**HPCPS A4392**

| Stoma Size | Plateau Size | Opaque w/ Easy-View Window |
|--------------|--------------|----------------------------|
| 9/16" (15mm) | V3 | 421629 |
| 9/16" (15mm) | V2 | 421627 |
| 9/16" (15mm) | V1 | 421624 |
| 4/5" (20mm) | V1 | 421625 |

**ESTEEM™ + SOFT CONVEX POUCH****Urostomy | Cut-to-Fit | Convex**

- Esteem™ + Soft Convex Adhesive (Standard Wear)
- Without Tape Collar
- Convex Wafer Depth = 4mm
- Soft Tap Closure
- Pouch Length = 7 2/3"
- Includes Belt Tabs
- 2 Esteem™ + Soft Convex Night Drainage Adapters per box
- Quantity = 10

HPCPS A4392

| Cutting Range | Plateau Size | Opaque w/ Easy-View Window |
|------------------------|--------------|----------------------------|
| 3/8 - 1 1/8" (10-28mm) | V3 | 422562 |
| 3/5 - 1 3/5" (15-40mm) | V2 | 422557 |
| 3/4 - 1 7/8" (20-47mm) | V1 | 422551 |

Urostomy | Pre-Cut | Convex**HPCPS A4392**

| Stoma Size | Plateau Size | Opaque w/ Easy-View Window |
|---------------|--------------|----------------------------|
| 4/5" (20mm) | V3 | 422564 |
| 1" (25mm) | V3 | 422565 |
| 4/5" (20mm) | V2 | 422558 |
| 1" (25mm) | V2 | 422559 |
| 1 1/6" (30mm) | V2 | 422560 |
| 1" (25mm) | V1 | 422552 |
| 1 1/6" (30mm) | V1 | 422553 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS**ONE-PIECE SYSTEMS**

- Esteem Body™ Pouches
- Esteem™ + Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

ActiveLife® Pouches

ACTIVELIFE® POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- 1 Tail Closure per box
- Quantity = See table below

HCPCS A5061 | No Filter

| Cutting Range | Transparent Qty = 1 | Transparent Qty = 10 | Transparent Qty = 20 |
|------------------------|---------------------|----------------------|----------------------|
| 3/4 - 2 1/2" (19-64mm) | 025433 | 022771 | 125330 |



ACTIVELIFE® POUCH

Cut-to-Fit | Drainable | Clipped Closure

- Durahesive® Plus (Extended Wear)
- Acrylic Tape Collar
- Pouch Length = 12"
- 1 Tail Closure per box
- Quantity = 10

HCPCS A4388 | No Filter

| Cutting Range | Transparent | Opaque |
|------------------------|-------------|--------|
| 3/4 - 2 1/2" (19-64mm) | 400598 | 400599 |



ACTIVELIFE® POUCH

Pre-Cut | Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Acrylic Tape Collar
- Pouch Length = 12" (*pouch length = 10")
- 1 Tail Closure per box
- Quantity = See table below



HCPCS A5061 | No Filter

| Stoma Size | Transparent Qty = 10 | Opaque Qty = 10 | Opaque Qty = 20 | Opaque* Qty = 10 | Opaque* Qty = 20 |
|---------------|----------------------|-----------------|-----------------|------------------|------------------|
| 3/4" (19mm) | 022764 | 022757 | 125338 | 022750* | 125331* |
| 1" (25mm) | 022765 | 022758 | 125339 | 022751* | N/A |
| 1 1/4" (32mm) | 022766 | 022759 | 125340 | 022752* | 125333* |
| 1 1/2" (38mm) | 022767 | 022760 | 125341 | 022753* | 125334* |
| 1 3/4" (45mm) | 022768 | 022761 | 125342 | 022754* | 125335* |
| 2" (50mm) | 022769 | 022762 | 125343 | 022755* | 125336* |
| 2 1/2" (64mm) | 022770 | 022763 | 125344 | 022756* | 125337* |

ACTIVELIFE® POUCH

Pre-Cut | Drainable | Convex | Clipped Closure

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Convex Wafer Depth = 5mm
- Pouch Length = 12"
- Includes Belt Tabs
- 1 Tail Closure per box
- Quantity = See table below



HCPCS A4390 | No Filter

| Stoma Size | Transparent Qty = 5 |
|---------------|---------------------|
| 3/4" (19mm) | 175777 |
| 7/8" (22mm) | 175778 |
| 1" (25mm) | 175779 |
| 1 1/8" (28mm) | 175780 |
| 1 1/4" (32mm) | 175781 |
| 1 3/8" (35mm) | 175782 |
| 1 1/2" (38mm) | 175783 |
| 1 3/4" (45mm) | 175784 |
| 2" (50mm) | 175785 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ONE-PIECE SYSTEMS

- Esteem Body™ Pouches
- Esteem™+ Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

ACTIVELIFE® POUCH

Cut-to-Fit | Stoma Cap

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Quantity = 30



HCPCS A5055 | Filter (Filter covers included)

| Cutting Range | Opaque |
|----------------------|--------|
| 3/4 - 2" (19 - 50mm) | 175611 |

ACTIVELIFE® POUCH

Urostomy | Pre-Cut | Accuseal™ Tap with Valve

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Anti-Reflux Valve
- Pouch Length = 8"
- Includes Belt Tabs
- 1 Accuseal Night Drainage Adapter per box
- Quantity = 10



HCPCS A4428

| Stoma Size | Transparent |
|---------------|-------------|
| 3/4" (19mm) | 650828 |
| 7/8" (22mm) | 650829 |
| 1" (25mm) | 650830 |
| 1 1/8" (28mm) | 650831 |
| 1 1/4" (32mm) | 650832 |
| 1 1/2" (38mm) | 650833 |

ACTIVELIFE® POUCH

Urostomy | Cut-to-Fit | Fold Over Tap

- Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 9"
- Includes Belt Tabs
- 2 Universal Night Drainage Adapters per box
- Quantity = 10



HCPCS A5071

| Cutting Range | Transparent |
|------------------------|-------------|
| 3/4 - 1 3/4" (19-45mm) | 064927 |

ACTIVELIFE® POUCH

Urostomy | Pre-Cut | Convex | Fold Over Tap

- Durahesive® (Extended Wear)
- Acrylic Tape Collar
- Convex Wafer Depth = 5mm
- Pouch Length = 8"
- Includes Belt Tabs
- 1 Universal Night Drainage Adapter per box
- Quantity = See table below



HCPCS A4393

| Stoma Size | Transparent Qty = 5 | Transparent Qty = 10 |
|---------------|---------------------|----------------------|
| 1/2" (13mm) | 175790 | 125361 |
| 5/8" (16mm) | 175791 | 125362 |
| 3/4" (19mm) | 175792 | 125363 |
| 7/8" (22mm) | 175793 | 125364 |
| 1" (25mm) | 175794 | 125365 |
| 1 1/8" (28mm) | 175795 | 125366 |
| 1 1/4" (32mm) | 175796 | 125367 |
| 1 3/8" (35mm) | 175797 | 125368 |
| 1 1/2" (38mm) | 175798 | 125369 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ONE-PIECE SYSTEMS

- Esteem Body™ Pouches
- Esteem™+ Pouches
- ActiveLife® Pouches

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

ACUTE KITS



Convatec's Acute Kits provide simple post-operative solutions designed to get your patient into the right system after surgery.

All of Convatec's Kits include NATURA™ Skin Barriers and NATURA™ Pouches and feature Durahesive® to provide initial tack, strong adhesion and extended wear time.

Systems include Cut-to-Fit and Moldable barrier options, as well as flat and low-pressure Accordion flanges.

Color Key & Sizing Guide







| COLOR MATCH |  |  |  |  |  |  |
|-------------|---|---|---|---|---|---|
| FLANGE SIZE | 1 1/4" (32mm) | 1 1/2" (38mm) | 1 3/4" (45mm) | 2 1/4" (57mm) | 2 3/4" (70mm) | 4" (100mm) |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

ACUTE KITS

- Drainable Flat
- Drainable Low-Pressure Accordion
- Urostomy Low-Pressure Accordion

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Drainable Flat

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Cut-to-Fit | Drainable | Invisiclose™ Tail Closure

- Durahesive® (Extended Wear)
- Without Tape Collar
- Pouch Color: Transparent
- Quantity = 5 Kits

HPCPS N/A | Non-Sterile

| | Cutting Range (mm) | Kit Product Code | KIT INCLUDES: | |
|---|--------------------|------------------|---------------|--------|
| | | | Skin Barrier | Pouch |
| ● | 13-45 | 423556 | 413156 | 416421 |
| ● | 13-57 | 423555 | 413157 | 416424 |
| ● | 13-89 | 423554 | 401905 | 416472 |

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Moldable | Drainable | InvisiClose™ Tail Closure

- Durahesive® (Extended Wear)
- Hydrocolloid Tape Collar (White)
- Pouch Color: Transparent
- Quantity = 5 Kits

HPCPS N/A | Non-Sterile

| | Stoma Size (mm) | Kit Product Code | KIT INCLUDES: | |
|---|-----------------|------------------|---------------|--------|
| | | | Skin Barrier | Pouch |
| ● | 22-33 | 432559 | 411802 | 416418 |
| ● | 33-45 | 432558 | 411804 | 416421 |
| ● | 45-56 | 423557 | 411806 | 416424 |

Drainable & Urostomy Low-Pressure Accordion

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Cut-to-Fit | Accordion Flange | Drainable | Invisiclose™ Tail Closure

- Durahesive® (Extended Wear)
- Acrylic Tape Collar (White)
- Pouch Color: Transparent
- Quantity = 5 Kits

HPCPS N/A | Non-Sterile

| | Cutting Range (mm) | Kit Product Code | KIT INCLUDES: | |
|---|--------------------|------------------|---------------|--------|
| | | | Skin Barrier | Pouch |
| ● | 13-45 | 423548 | 421462 | 416424 |

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

Moldable | Accordion Flange | Drainable | Invisiclose™ Tail Closure

- Durahesive® (Extended Wear)
- Hydrocolloid Tape Collar (Tan)
- Pouch Color: Transparent
- Quantity = 5 Kits

HPCPS N/A | Non-Sterile

| | Stoma Size (mm) | Kit Product Code | KIT INCLUDES: | |
|---|-----------------|------------------|---------------|--------|
| | | | Skin Barrier | Pouch |
| ● | 33-45 | 423549 | 421041 | 416424 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

ACUTE KITS

- Drainable Flat
- Drainable Low-Pressure Accordion
- Urostomy Low-Pressure Accordion

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

**Cut-to-Fit | Accordion | Urostomy |
Accuseal™ Tap with Valve**

- Available in Stomahesive® (Standard Wear) and Durahesive® (Extended Wear)
- Acrylic Tape Collar (White)
- 1 Accuseal Night Drainage Adapter per box
- Pouch Color: Transparent
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

| | Cutting Range (mm) | Kit Product Code | KIT INCLUDES: | |
|---|--------------------|------------------|---------------|--------|
| | | | Skin Barrier | Pouch |
| ● | 13-33 | 423551 | 421460 | 401545 |
| ● | 13-45 | 423550 | 421462 | 401546 |

NATURA™ SKIN BARRIER & NATURA™ + POUCH KIT

**Moldable | Accordion | Urostomy |
Accuseal™ Tap with Valve**

- Durahesive® (Extended Wear)
- Hydrocolloid Tape Collar (Tan)
- Anti-Reflux Valve
- 1 Accuseal Night Drainage Adapter per box
- Pouch Color: Transparent
- Quantity = 5 Kits

HCPCS N/A | Non-Sterile

| | Stoma Size (mm) | Kit Product Code | KIT INCLUDES: | |
|---|-----------------|------------------|---------------|--------|
| | | | Skin Barrier | Pouch |
| ● | 22-33 | 423553 | 421040 | 401545 |
| ● | 33-45 | 423552 | 421041 | 401546 |

CONVATEC LOOP OSTOMY ROD

Sterile

- Quantity = 10

HCPCS A4421

| | Flange Size | Size | Product Code |
|---|---------------|------|--------------|
| ● | 2 3/4" (70mm) | 65mm | 022355 |
| ● | 4" (100mm) | 90mm | 022356 |



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

ACUTE KITS

- Drainable Flat
- Drainable Low-Pressure Accordion
- **Urostomy Low-Pressure Accordion**

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

PEDIATRIC





Products



Little Ones® is a pediatric Two-piece and One-piece system designed for babies smaller than 10 lbs, and children over 10 lbs.

The pouches have soft-welded edges and smaller pouch sizes designed to provide extra comfort for little ones.

Color Key & Sizing Guide

| COLOR MATCH |  |  |  |  |
|-----------------------|---|---|---|---|
| FLANGE/ COUPLING SIZE | 0 - 9/10" (0-23mm) | 1/5 - 1 1/4" (4.8-31mm) | 1 1/4" (32mm)* | 1 3/4" (45mm)* |

*Indicates Flange size measurement



When ordering a Pediatric Two-Piece System, **make sure you select the same flange/ coupling size for the skin barrier and pouch.** Use the colored dots in the table to help.

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

PEDIATRIC PRODUCTS

- Extra-Small Systems
- Two-Piece Systems
- One-Piece Systems

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

Little Ones® Extra Small Two-Piece Ostomy System

For babies smaller than 10lbs/4.5kg


LITTLE ONES® EXTRA SMALL SKIN BARRIER

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Adhesive Coupling Technology™
- No Starter Hole
- Quantity = 5

HCPCS A4414

| Cutting Range | Without Tape Collar |
|----------------------|---------------------|
| ● 0 - 9/10" (0-23mm) | 411632 |



LITTLE ONES® EXTRA SMALL POUCH

Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Adhesive Coupling Technology™
- Pouch Length = 6"
- Quantity = 10

HCPCS A5063 | No Filter

| Stoma Size | Transparent |
|----------------------|-------------|
| ● 0 - 9/10" (0-23mm) | 411631 |



Little Ones® Extra Small One-Piece Ostomy System

For babies smaller than 10lbs/4.5kg

LITTLE ONES® EXTRA SMALL POUCH

Cut-to-Fit | Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Without Tape Collar
- No Starter Hole
- Pouch Length = 6"
- Quantity = 10

HCPCS A5061 | No Filter

| Stoma Size | Transparent |
|--------------------|-------------|
| 0 - 9/10" (0-23mm) | 411633 |



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC
PRODUCTS

PEDIATRIC PRODUCTS

- Extra-Small Systems
- Two-Piece Systems
- One-Piece Systems

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample


Little Ones® Two-Piece Ostomy System

Suggested for children over 10lbs/4.5kg

LITTLE ONES® SKIN BARRIER

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Adhesive Coupling Technology™
- Quantity = 5



HCPCS A4414

| | Cutting Range | Hydrocolloid Tape Collar WHITE | No Tape Collar |
|---|-------------------------|-----------------------------------|----------------|
| ● | 1/5 - 1 1/4" (4.8-31mm) | 411642 | 411641 |

LITTLE ONES® POUCH

Drainable | InvisiClose™ Tail Closure

- Stomahesive® (Standard Wear)
- Adhesive Coupling Technology™
- Pouch Length = 8"
- Quantity = 10



HCPCS A4425 | Filter (Filter covers included)

| | Stoma Size | Transparent | Opaque |
|---|--------------------------|-------------|--------|
| ● | 1/5 - 1 1/4" (4.8-31 mm) | 411637 | 411638 |

LITTLE ONES® POUCH

Closed-End

- Stomahesive® (Standard Wear)
- Adhesive Coupling Technology™
- Pouch Length = 6"
- Quantity = 10



HCPCS A4419 | Filter (Filter covers included)

| | Stoma Size | Transparent | Opaque |
|---|--------------------------|-------------|--------|
| ● | 1/5 - 1 1/4" (4.8-31 mm) | 411640 | 411639 |

LITTLE ONES® SKIN BARRIER

Cut-to-Fit

- Stomahesive® (Standard Wear)
- Flange System
- Quantity = 5



HCPCS A4414

| | Flange Size | Cutting Range | Acrylic Tape Collar WHITE |
|---|---------------|--------------------------|------------------------------|
| ● | 1 1/4" (32mm) | 1/2 - 3/4" (13-19mm) | 401925 |
| ● | 1 3/4" (45mm) | 1 1/2 - 1 1/4" (13-32mm) | 401926 |

LITTLE ONES® POUCH

Drainable | Clipped Closure

- Stomahesive® (Standard Wear)
- Flange System
- Pouch Length = 6"
- Includes Belt Tabs
- 10 Soft, Moldable Tail Closures per box (Designed for Pediatric Use)
- Quantity = 10



HCPCS A5063 | No Filter

| | Flange Size | Transparent |
|---|---------------|-------------|
| ● | 1 1/4" (32mm) | 401927 |
| ● | 1 3/4" (45mm) | 401928 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC
PRODUCTS

PEDIATRIC PRODUCTS

- Extra-Small Systems
- **Two-Piece Systems**
- One-Piece Systems

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

LITTLE ONES® POUCH

Closed-End

- Flange System
- Pouch Length = 5"
- Includes Belt Tabs
- Quantity = 20

**HCPCS A5054 | No Filter**

| | Flange Size | Opaque |
|---|---------------|--------|
| ● | 1 1/4" (32mm) | 401930 |
| ● | 1 3/4" (45mm) | 401931 |

LITTLE ONES® POUCH

Urostomy | Fold Over Tap

- Flange System
- Pouch Length = 6"
- Includes Belt Tabs
- 1 Universal Night Drainage Adapter per box
- Quantity = 10

**HCPCS A5073 | No Filter**

| | Flange Size | Transparent |
|---|---------------|-------------|
| ● | 1 1/4" (32mm) | 401929 |

Little Ones® One-Piece Ostomy System

Suggested for children over 10lbs/4.5kg

LITTLE ONES® POUCH

**Cut-to-Fit | Drainable |
Clipped Closure**

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 6"
- 15 Soft, Moldable Tail Closures per box (Designed for Pediatric Use)
- Quantity = 15

**HCPCS A5061 | No Filter**

| Cutting Range | Transparent |
|--------------------|-------------|
| 5/16 - 2" (8-50mm) | 020922 |

LITTLE ONES® POUCH

**Urostomy | Cut-to-Fit |
Fold Over Tap**

- Modified Stomahesive® (Standard Wear)
- Without Tape Collar
- Pouch Length = 5"
- 2 Universal Night Drainage Adapters per box
- Quantity = 15

**HCPCS A5071**

| Cutting Range | Transparent |
|--------------------|-------------|
| 5/16 - 1" (8-25mm) | 020917 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC
PRODUCTS**PEDIATRIC PRODUCTS**

- Extra-Small Systems
- Two-Piece Systems
- **One-Piece Systems**

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

ACCESSORIES

CHOOSING THE RIGHT ACCESSORIES

Whether you are looking for better skin protection or added security to your ostomy system, with the right combination of accessories, you can feel even more comfortable and at ease.



Skin Preparation



Pouching Support



Added Security



Pouch Removal

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

ACCESSORIES

- Skin Preparation
- Pouching Support
- Added Security
- Pouch Removal
- Irrigation

WOUND & FISTULA

RESOURCES

Request Product Sample



Skin Preparation

Solutions that help promote healthy skin by reducing irritation and damage.

ESENTA™

Sting-Free Skin Barrier Wipe

- Creates a breathable barrier on skin
- Protects skin for up to 72 hours¹⁸
- Reduces the risk of skin damage and irritation caused by body wastes, enzymes, and adhesives
- Safe and effective even with frequent use
- Silicone based
- Alcohol free
- CHG Compatible¹⁹



HCPCS A5120

| Liquid per Wipe | Product Code | Box of | Boxes per Case* |
|-----------------|--------------|--------|-----------------|
| 3 mL | 423392 | 25 | 20 |

*Sold by the case

ESENTA™

Sting-Free Skin Barrier Pump Spray

HCPCS A4369

| Spray Bottle Size | Product Code | Box of | Boxes per Case* |
|-------------------|--------------|--------|-----------------|
| 28 mL | 423286 | 1 | 12 |
| 50 mL | 423288 | 1 | 12 |

*Sold by the case

Continued ➤

ESENTA™

Sting-Free Skin Barrier Foam Applicator



HCPCS A5120

| Liquid per Applicator | Product Code | Box of | Boxes per Case* |
|-----------------------|--------------|--------|-----------------|
| 1 mL | 423284 | 25 | 4 |
| 1 mL | 423283 | 5 | 20 |
| 3 mL | 423285 | 25 | 4 |
| 3 mL | 423287 | 5 | 20 |

SION BIOTEXT®

Protective Barrier Wipe



- Designed to create a protective layer that shields skin from damage
- Helps reduce skin stripping and irritation from body waste
- May help increase the contact or adhesion to skin
- Contains Alcohol

HCPCS A5120

| Product Code | Box of |
|--------------|--------|
| 423779 | 50 |
| 423784 | 100 |

STOMAHESIVE®

Protective Powder

- Helps protect skin by absorbing excess moisture
- For use on moist skin around a stoma
- Lightly dust affected area with powder, wiping off excess. May be followed up with a skin barrier wipe or spray



HCPCS A4371

| Size | Product Code | Box of |
|-------|--------------|----------|
| 1 oz. | 025510 | 1 bottle |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

ACCESSORIES

- Skin Preparation
- Pouching Support
- Added Security
- Pouch Removal
- Irrigation

WOUND & FISTULA

RESOURCES

Request Product Sample

2 Pouching Support

Products that help create the best fit around the stoma, prevent leakage and help level out the application surface area for better pouch performance.

EAKIN COHESIVE®

Seals

- Moldable ring that helps seal the area around hard-to-fit stomas from any leakage under the skin barrier
- Acts as a filler on uneven skin surfaces to help create a flatter area for the skin barrier to adhere to
- Designed to absorb moisture and prevent it coming in contact with, and causing harm to the skin—forms a gel that protects the skin

HCPCS A4385

| Size | Product Code | Depth | Box of |
|------------------|--------------|--------|--------|
| Small 2" | 839002 | 4.2 mm | 20 |
| Slims 2" | 839005 | 3 mm | 10 |
| Large 4" | 839001 | 2.1 mm | 10 |
| StomaWrap 3 1/3" | 839006 | 3 mm | 10 |



WATCH NOW

to learn more about
eakin Cohesive® Seals



STOMAHESIVE®

Paste

- Helps to protect exposed skin between the base of the stoma and the opening of the skin barrier. Can also be used to fill crevices, skin folds and scars to help create a flat pouching surface for optimal appliance attachment
- Helps to prevent leakage under the skin barrier that may result in skin breakdown
- Pectin-based
- Contains alcohol

HCPCS A4406

| Size | Product Code | Box of |
|-------|--------------|--------|
| 2 oz. | 183910 | 1 tube |



EAKIN COHESIVE®

Paste

- Use to fill crevices, skin folds and scars to help create an even pouching surface
- Leakage prevention by helping form a waterproof barrier around the stoma
- Alcohol free

HCPCS A4405

| Size | Product Code | Box of |
|---------|--------------|--------|
| 2.1 oz. | 839010 | 1 tube |



TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

ACCESSORIES

- Skin Preparation
- **Pouching Support**
- Added Security
- Pouch Removal
- Irrigation

WOUND & FISTULA

RESOURCES

Request Product Sample

STOMAHESIVE®

Skin Barrier

- Helps fill uneven skin surfaces to create a flat pouching surface for optimal appliance attachment
- Designed to help protect the skin from contact with bodily fluids
- Large cutting surface to accommodate large and swollen stomas

HCPCS A4362

| Size | Product Code | Box of |
|--------|--------------|--------|
| 4 x 4" | 021712 | 5 |

HCPCS A5122

| Size | Product Code | Box of |
|--------|--------------|--------|
| 8 x 8" | 021715 | 3 |

EAKIN COHESIVE®

Skin Barrier

- Designed to help protect the skin from contact with body fluids
- Acts as a filler on uneven skin surfaces to help create a flat pouching surface for optimal appliance attachment
- Moldable and moisture absorbing

HCPCS A4385

| Size | Product Code | Box of |
|--------|--------------|--------|
| 4 x 4" | 839004 | 5 |
| 4 x 8" | 839003 | 5 |

STOMAHESIVE®

Strips

- Helps to create a secure seal between the stoma and the appliance, to protect the skin against stoma output
- Can be cut, torn and shaped into a variety of shapes to help prevent leakage around the stoma by acting as a filler
- Designed to absorb moisture from the stoma, helping to create a snug fit around the stoma

HCPCS A4362

| Size | Product Code | Box of |
|--------------------------|--------------|-----------|
| 3 cards of 5 strips each | 025542 | 15 strips |

SUR-FIT® NATURA™

Disposable Convex Inserts

- Convex inserts that can be added to a flat skin barrier to create a shallow convexity—insert into the inner diameter of the plastic flange on the pouch facing side of a SUR-FIT® Natura™ Durahesive® and Stomahesive® skin barrier
- Select an insert that is at least $\frac{1}{4}$ " larger than the maximum diameter of the stoma
- For SUR-FIT® Natura™ skin barriers only. Not for use with Convatec Moldable Technology™ skin barriers, or any other product
- Quantity = 5

HCPCS A5093

| | For Use With Flange Size | Stoma Size no Wider Than | Convex Insert Opening Size | Product Code |
|---|--------------------------|--------------------------|----------------------------|--------------|
| ● | 1 1/2" (38mm) | 1/2" | 3/4" (19mm) | 404006 |
| ● | 1 1/2" (38mm) | 3/4" | 1" (25mm) | 404008 |
| ● | 1 3/4" (45mm) | 7/8" | 1 1/8" (29mm) | 404009 |
| ● | 1 3/4" (45mm) | 1" | 1 1/4" (32mm) | 404010 |
| ● | 1 3/4" (45mm) | 1 1/8" | 1 3/8" (35mm) | 404011 |
| ● | 2 1/4" (57mm) | 1 1/4" | 1 1/2" (38mm) | 404012 |
| ● | 2 1/4" (57mm) | 1 3/8" | 1 5/8" (41mm) | 404013 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA

RESOURCES

ACCESSORIES

- Skin Preparation
- **Pouching Support**
- Added Security
- Pouch Removal
- Irrigation

Request Product Sample



Added Security

An extra layer of protection to keep your ostomy system in place.

ESENTA™

Lubricating Deodorant

- Lubricating spray coats the pouch to prevent stool build up around the stoma that can compromise the seal
- Stops odor before it starts
- Made with natural Eucalyptus, lavender, pine oils and other natural ingredients
- Fast, easy spray covers the pouch thoroughly without the mess



HCPCS A4394

| Spray Bottle Size | Product Code | Box of | Boxes per Case* |
|-------------------|--------------|--------|-----------------|
| 236mL | 423574 | 1 | 25 |

*Sold by the case

EASE™

Strips

- Helps keep your skin barrier in place
- Can be applied to cover the edges of the skin barrier
- Designed to be thin and flexible
- Water-resistant



HCPCS A4362 (4 strips = 1 billable unit of A4362)

| Size | Product Code | Quantity |
|--------------------------|--------------|-----------|
| 1 1/5 x 3 1/2" (3 x 9cm) | 422163 | 80 strips |

CONVATEC®

Ostomy Appliance Belt

- Used to help keep your ostomy system in place



HCPCS A4367

| Size | Product Code | Box of |
|----------------------|--------------|--------|
| Adjustable up to 42" | 175507 | 1 |

DIAMONDS™

Gelling and Odor Control Sachets

- Gels liquid contents
- Contains activated charcoal which helps to eliminate odors



HCPCS A4422

| Product Code | Box of |
|--------------|-------------|
| TR105 | 100 Sachets |

DUOLOCK™

Curved Tail Closure

- Helps secure the end of Convatec® Drainable pouch



HCPCS A4363

| Product Code | Box of |
|--------------|--------|
| 175652 | 10 |

CONVATEC®

Tail Closure

- Helps secure the end of Convatec® Drainable pouch



HCPCS A4363

| Product Code | Box of |
|--------------|--------|
| 175635 | 10 |

CONVATEC®

Night Drainage Container Set*

- Plastic vessel to use for overnight urine collection
- Container capacity of 2000 cc/ml
- Reusable, can be washed out between uses



HCPCS A5102

| Product Code | Quantity |
|--------------|----------|
| 027060 | 1 set |

*Contains: 1 Convatec® Night Drainage Container with tubing.
1 Convatec® Night Drainage Container Cover and 1 Universal Adapter

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA

RESOURCES

ACCESSORIES

- Skin Preparation
- Pouching Support
- **Added Security**
- Pouch Removal
- Irrigation

Request Product Sample

Added Security Cont'd

CONVATEC®



Night Drainage Container Tubing

- Replacement tubing for the Convatec® Night Drainage Container

HCPCS A4331

| Product Code | Box of |
|--------------|--------|
| 027062 | 1 |



Pouch Removal

Solutions that help promote healthy skin by protecting it from irritation and damage caused by adhesives.

ESENTA™



Sting-Free Adhesive Remover Wipe

- For fast and easy adhesive release
- Reduces painful skin stripping and skin damage
- Quickly removes adhesive residue
- Feels silky and smooth, no residue left behind
- Silicone based
- Alcohol free
- CHG Compatible¹⁹

HCPCS A4456

| Liquid per Wipe | Product Code | Box of | Boxes per Case* |
|-----------------|--------------|--------|-----------------|
| 3 mL | 423391 | 25 | 20 |

*Sold by the case

ESENTA™

Sting-Free Adhesive Remover Spray

- For fast and easy adhesive release
- Reduces painful skin stripping and skin damage
- Quickly removes adhesive residue
- Feels silky and smooth, no residue left behind
- Silicone based
- Alcohol free
- CHG Compatible¹⁹



HCPCS A4455

| Spray Bottle Size | Product Code | Box of | Boxes per Case* |
|-------------------|--------------|--------|-----------------|
| 50 mL | 423289 | 1 | 12 |
| 150 mL | 423290 | 1 | 12 |

*Sold by the case

SION BIOTEXT®



Adhesive Remover Wipe

- Helps gently remove adhesives, including skin barriers, tapes and hydrocolloid dressings
- Contains alcohol

HCPCS A4456

| Product Code | Box of |
|--------------|--------|
| 423780 | 50 |
| 423783 | 100 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA

RESOURCES

ACCESSORIES

- Skin Preparation
- Pouching Support
- Added Security
- **Pouch Removal**
- Irrigation

Request Product Sample

Irrigation

VISI-FLOW®

Irrigator with Stoma Cone

HCPCS A4398 & A4399

Product Code

Box of

401989

1



VISI-FLOW®

Irrigation Sleeve

- Compatible with Natura™ and SUR-FIT® Natura™ skin barriers
- Color: Transparent
- 1 Tail Closure per box
- Quantity = 5



HCPCS A4436

Flange Size

Product Code

| | | |
|---|---------------|--------|
| ● | 1 1/2" (38mm) | 401911 |
| ● | 1 3/4" (45mm) | 401912 |
| ● | 2 1/4" (57mm) | 401913 |
| ● | 2 3/4" (70mm) | 401914 |

VISI-FLOW®

Stoma Cone

HCPCS A4399

Product Code

Box of

022736

1



VISI-FLOW®

Irrigation Sleeve Tail Closure

HCPCS A4363

Product Code

Box of

175650

10



VISI-FLOW®

Irrigation Adaptor Faceplate

- Compatible with Visi-Flow® Irrigation Sleeves
- Quantity = 1



HCPCS A4361

Flange Size

Product Code

| | | |
|---|---------------|--------|
| ● | 1 3/4" (45mm) | 401918 |
| ● | 2 3/4" (70mm) | 401919 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA

RESOURCES

ACCESSORIES

- Skin Preparation
- Pouching Support
- Added Security
- Pouch Removal
- Irrigation

Request Product Sample

eakin® WOUND & FISTULA Pouches

eakin® wound pouches can help contain fistula output and handle exuding wounds. Pouch options include the ability to connect tubing for remote drainage, or a fold & tuck option.



WATCH NOW

to learn how to
apply an eakin®
Wound Pouch

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

EAKIN®

Wound Pouch (1.2 x 1.8" Cutting Area)

- Quantity = 10

HCPCS A6154

| Cutting Range | Fold and Tuck Closure |
|------------------------|-----------------------|
| 1.2 x 1.8" (3 x 4.5cm) | 839250 |



EAKIN®

Wound Pouch (3.14 X 2.4" Cutting Area)

- Quantity = 10

HCPCS A6154

| Cutting Range | Remote Drainage Attachment and Tap Closure |
|-------------------------|--|
| 3.14 x 2.4" (8 x 6.1cm) | 839267 |



EAKIN®

Wound Pouch (3.0 X 4.3" Cutting Area)

- Quantity = 10

HCPCS A6154

| Cutting Range | Remote Drainage Attachment and Tap Closure | Fold and Tuck Closure |
|---------------------------|--|-----------------------|
| 3.0 x 4.3" (7.6 x 10.9cm) | 839261 | 839251 |



EAKIN®

Wound Pouch (4.3 X 6.9" Cutting Area)

- Quantity = 10

HCPCS A6154

| Cutting Range | Remote Drainage Attachment and Tap Closure | Fold and Tuck Closure |
|----------------------------|--|-----------------------|
| 4.3 x 6.9" (10.9 x 17.5cm) | 839262 | 839252 |



EAKIN®

Horizontal Cutting Wound Pouch (9.7 X 6.3" Cutting Area)

- Quantity = 5

HCPCS A6154

| Cutting Range | Remote Drainage Attachment and Tap Closure with Access Window | Fold and Tuck Closure |
|--------------------------|---|-----------------------|
| 9.7 x 6.3" (24.6 x 16cm) | 839264 | 839253 |



EAKIN®

Vertical Cutting Wound Pouch (6.3 X 9.7" Cutting Area)

- Quantity = 5

HCPCS A6154

| Cutting Range | Remote Drainage Attachment and Tap Closure | Remote Drainage Attachment and Tap Closure with Access Window |
|--------------------------|--|---|
| 6.3 x 9.7" (16 x 24.6cm) | 839265 | 839268 |



EAKIN®

Vertical Cutting Wound Pouch (5.1 X 11.4" Cutting Area)

- Quantity = 5

HCPCS A6154

| Cutting Range | Remote Drainage Attachment and Tap Closure |
|-------------------------|--|
| 5.1 x 11.4" (13 x 29cm) | 839266 |



EAKIN®

Eakin Access Window

- Add to Eakin wound and fistula pouches to allow instant access to the wound (may be used with 839251, 839252, 839253, 839262, 839265, 839266)
- Quantity = 5

HCPCS n/a

| Size | Remote Drainage Attachment and Tap Closure |
|------|--|
| N/A | 839280 |



me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample

RESOURCES

Our resources will help you support your patients living with an ostomy.



- ✓ Patient Ostomy Support & Resources
- ✓ QR Code Quick Guide
- ✓ Medicare Utilization Charts
- ✓ Product Code Index
- ✓ References

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FIST

RESOURCES

RESOURCES

- Patient Ostomy Support & Resources
- Quick Links Guide
- Medicare Utilization Chart
- Product Code Index
- References

Request Product Sample

OSTOMY RESOURCES & SUPPORT

for You and Your Patient

Convatec

Convatec is a global medical products and technologies company developing ostomy products and solutions for almost 40 years.

Toll Free: 1-800-422-8811
Convatec.com/ostomy

Ostomysecrets®

A Convatec company, offers a full line of stylish and functional undergarments, swim wear and accessories designed to support your ostomy solution.

Toll Free: 1-877-613-6246
ostomysecrets.com

me+™

Support, insights & product information for those living with an ostomy. Enrolling in me+™ is fast, easy and free.

Toll Free: 1-800-422-8811
e-mail: CIC@Convatec.com
mepluscare.com

Ostomy Nurse Solutions™ App

A fast and easy way to help you support your patient without leaving the bedside.

ostomynursesolutions.com

My Ostomy Journey™ App

Patients can conveniently track ostomy-related activities without the hassle of keeping up with traditional paper check lists or diaries

myostomyjourney.com

OTHER OSTOMY RESOURCES

Wound Ostomy Continence Nurses Society™ (WOCN®)

Professional nursing society dedicated to providing expert health care to individuals with wound, ostomy and continence concerns.

Toll Free: 1-888-224-9626
wocn.org

American Cancer Society

Nationwide community-based, voluntary healthcare organization.

Toll Free: 1-800-ACS-2345
cancer.org

The United Ostomy Associations of America (UOAA)

National organization providing support, information and advocacy.

Toll Free: 1-800-826-0826
ostomy.org

Crohn's & Colitis Foundation of America

Non-profit, volunteer-driven organization dedicated to finding a cure for Crohn's Disease and ulcerative colitis.

Toll Free: 1-800-932-2423
ccfa.org

OSTOMY RESOURCES & SUPPORT

Ostomy Nurse Solutions™ App

My Ostomy Journey™ App

Application Videos

Convatec Moldable Technology™
Product Information

me+™ Patient Enrollment

Accessories

Convexity Product Information

Discontinued Product Information

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FIST

RESOURCES

RESOURCES

- Patient Ostomy Support & Resources
- Quick Links Guide
- Medicare Utilization Chart
- Product Code Index
- References

Request Product Sample

| HCPCS | HCPCS Descriptors | Medicare Usual Maximum Qty. |
|-------|---|-----------------------------|
| A4331 | Extension Drainage Tube | 1 per 6 months |
| A4361 | Ostomy Faceplate | 3 per 6 months |
| A4362 | Skin Barrier, Solid, 4x4 or Equivalent | 20 per month |
| A4363 | Ostomy Clamp, Any Type, Replacement Only | As needed |
| A4367 | Ostomy Belt | 1 per month |
| A4369 | Ostomy Skin Barrier Liquid (spray or brush) | 2 oz per month |
| A4371 | Ostomy Skin Barrier Powder | 10 oz per 6 months |
| A4385 | Skin Barrier, Solid 4x4 or Equivalent, Extended Wear | As needed |
| A4388 | Pouch, Drainable, w/Extended Wear Barrier Attached | As needed |
| A4389 | Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) | 20 per month |
| A4390 | Pouch, Drainable, w/Extended Wear Barrier, Attached, w/convexity | As needed |
| A4392 | Ostomy Pouch, Urinary, w/ Barrier Attached, w/ Built-In Convexity (1-piece) | 20 per month |
| A4393 | Ostomy Pouch, Urinary, w/ Extended Wear Barrier Attached, w/ Built-In Convexity (1-piece) | 20 per month |
| A4394 | Ostomy Deodorant (with or without lubricant) | 8oz per month |
| A4436 | Irrigation Sleeve | 4 per month |
| A4398 | Ostomy Irrigation Bag | 2 per 6 months |
| A4399 | Ostomy Irrigation Cone/Catheter, and Brush | 2 per 6 months |
| A4402 | Lubricant | 4 oz per month |
| A4405 | Ostomy Skin Barrier, Non-Pectin Based, Paste | 4 oz per month |
| A4406 | Ostomy Skin Barrier, Pectin-Based, Paste | 4 oz per month |
| A4407 | Skin Barrier, w/Flange, Extended Wear w/ Convexity | As needed |
| A4409 | Skin Barrier, w/Flange, Extended Wear <= 4x4 | As needed |
| A4410 | Skin Barrier, w/Flange, Extended Wear > 4x4 | As needed |
| A4413 | Pouch, Drainable, High Output, for Barrier w/Flange | As needed |
| A4414 | Skin Barrier, w/Flange, <= 4x4 | 20 per month |
| A4415 | Skin Barrier, w/Flange, > 4x4 | 20 per month |

| HCPCS | HCPCS Descriptors | Medicare Usual Maximum Qty. |
|-------|---|-----------------------------|
| A4416 | Pouch, closed, w/Barrier Attached, w/Filter | 60 per month |
| A4417 | Ostomy Pouch, Closed, w/ Barrier Attached, w/ Built-In Convexity, w/ Filter (1-piece) | 60 per month |
| A4419 | Pouch, closed for use on Barrier w/Flange, w/Filter | 60 per month |
| A4421 | Ostomy Supply, Miscellaneous | As needed |
| A4422 | Ostomy Absorbent Material (Sheet/Pad/Crystal Packet) | No Fee Schedule |
| A4424 | Ostomy Pouch, Drainable, w/ Barrier Attached, w/ Filter (1-piece) | 20 per month |
| A4425 | Pouch, Drainable, with filter for use w/Flange | 20 per month |
| A4428 | Pouch, Urinary, w/ Extended Wear Barrier w/Faucet-type Tap | As needed |
| A4432 | Pouch, Urinary, for use on Barrier w/Flange, w/Faucet-type Tap | 20 per month |
| A4455 | Adhesive Remover or Solvent | 2 oz per month |
| A4456 | Adhesive Remover Wipes, Any Type | 50 per month |
| A5054 | Pouch Closed, for Barrier w/Flange | 60 per month |
| A5055 | Stoma Cap, each | 31 per month |
| A5056 | Pouch, Drainable, w/Extended Wear Barrier Attached, w/filter | 20 per month |
| A5057 | Pouch, Drainable, w/Extended Wear Barrier, w/convexity, w/filter | 20 per month |
| A5061 | Pouch, Drainable, w/Barrier Attached | 20 per month |
| A5063 | Pouch, Drainable, for use w/Flange | 20 per month |
| A5071 | Pouch, Urinary, w/Barrier Attached | 20 per month |
| A5073 | Pouch, Urinary, for Barrier w/Flange | 20 per month |
| A5093 | Ostomy Accessory, Convex Insert | 10 per month |
| A5102 | Bedside Drainage Bottle w/ or wo/ tubing, rigid or expandable | 2 per 6 months |
| A5120 | Skin Barrier, Wipes, each | 150 per 6 months |
| A5121 | Skin Barrier; Solid, 6 X 6 or Equivalent | 20 per month |
| A5122 | Skin Barrier, Solid, 8 X 8 or Equivalent | 20 per month |
| A6154 | Wound Pouch | As needed |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FIST

RESOURCES

RESOURCES

- Patient Ostomy Support & Resources
- Quick Links Guide
- **Medicare Utilization Chart**
- Product Code Index
- References

Request Product Sample

| SKU # | PG | SKU # | PG |
|--------|-----|--------|-----|
| 020917 | 88 | 125272 | 40 |
| 020922 | 88 | 125273 | 40 |
| 021712 | 95 | 125274 | 40 |
| 021715 | 95 | 125275 | 40 |
| 022355 | 80 | 125276 | 40 |
| 022356 | 80 | 125277 | 40 |
| 022736 | 101 | 125330 | 71 |
| 022750 | 72 | 125331 | 72 |
| 022751 | 72 | 125333 | 72 |
| 022752 | 72 | 125334 | 72 |
| 022753 | 72 | 125335 | 72 |
| 022754 | 72 | 125336 | 72 |
| 022755 | 72 | 125337 | 72 |
| 022756 | 72 | 125338 | 72 |
| 022757 | 72 | 125339 | 72 |
| 022758 | 72 | 125340 | 72 |
| 022759 | 72 | 125341 | 72 |
| 022760 | 72 | 125342 | 72 |
| 022761 | 72 | 125343 | 72 |
| 022762 | 72 | 125344 | 72 |
| 022763 | 72 | 125361 | 74 |
| 022764 | 72 | 125362 | 74 |
| 022765 | 72 | 125363 | 74 |
| 022766 | 72 | 125364 | 74 |
| 022767 | 72 | 125365 | 74 |
| 022768 | 72 | 125366 | 74 |
| 022769 | 72 | 125367 | 74 |
| 022770 | 72 | 125368 | 74 |
| 022771 | 71 | 125369 | 74 |
| 025433 | 71 | 175507 | 97 |
| 025510 | 92 | 175611 | 73 |
| 025542 | 96 | 175635 | 98 |
| 027060 | 98 | 175650 | 102 |
| 027062 | 99 | 175652 | 98 |
| 064927 | 74 | 175777 | 72 |
| 125257 | 39 | 175778 | 72 |
| 125258 | 39 | 175779 | 72 |
| 125259 | 39 | 175780 | 72 |
| 125260 | 39 | 175781 | 72 |
| 125261 | 39 | 175782 | 72 |
| 125262 | 39 | 175783 | 72 |
| 125263 | 39 | 175784 | 72 |
| 125264 | 39 | 175785 | 72 |
| 125265 | 39 | 175790 | 74 |
| 125266 | 39 | 175791 | 74 |
| 125267 | 40 | 175792 | 74 |
| 125268 | 40 | 175793 | 74 |
| 125269 | 40 | 175794 | 74 |
| 125270 | 40 | 175795 | 74 |
| 125271 | 40 | 175796 | 74 |

| SKU # | PG | SKU # | PG |
|--------|----|--------|-----|
| 175797 | 74 | 401576 | 39 |
| 175798 | 74 | 401577 | 39 |
| 183910 | 94 | 401905 | 35 |
| 400598 | 71 | 401906 | 39 |
| 400599 | 71 | 401909 | 46 |
| 401500 | 43 | 401911 | 101 |
| 401501 | 43 | 401912 | 101 |
| 401502 | 43 | 401913 | 101 |
| 401503 | 43 | 401914 | 101 |
| 401504 | 43 | 401918 | 102 |
| 401505 | 43 | 401919 | 102 |
| 401506 | 43 | 401925 | 86 |
| 401507 | 43 | 401926 | 86 |
| 401508 | 43 | 401927 | 86 |
| 401509 | 43 | 401928 | 86 |
| 401510 | 43 | 401929 | 87 |
| 401511 | 43 | 401930 | 87 |
| 401512 | 43 | 401931 | 87 |
| 401513 | 43 | 401932 | 43 |
| 401514 | 43 | 401933 | 43 |
| 401515 | 43 | 401934 | 43 |
| 401533 | 48 | 401935 | 43 |
| 401534 | 48 | 401936 | 43 |
| 401535 | 48 | 401989 | 101 |
| 401536 | 48 | 403948 | 49 |
| 401537 | 48 | 404006 | 96 |
| 401538 | 48 | 404008 | 96 |
| 401539 | 48 | 404009 | 96 |
| 401540 | 48 | 404010 | 96 |
| 401541 | 48 | 404011 | 96 |
| 401542 | 47 | 404012 | 96 |
| 401543 | 47 | 404013 | 96 |
| 401544 | 47 | 404014 | 43 |
| 401545 | 47 | 404015 | 43 |
| 401546 | 47 | 404016 | 43 |
| 401547 | 47 | 404017 | 43 |
| 401548 | 47 | 404018 | 43 |
| 401549 | 47 | 404025 | 43 |
| 401550 | 47 | 404026 | 43 |
| 401551 | 47 | 404027 | 43 |
| 401552 | 47 | 404028 | 43 |
| 401553 | 47 | 404029 | 43 |
| 401554 | 47 | 404030 | 43 |
| 401555 | 47 | 404031 | 43 |
| 401557 | 44 | 404032 | 43 |
| 401558 | 44 | 404033 | 43 |
| 401559 | 44 | 404034 | 43 |
| 401573 | 39 | 404592 | 38 |
| 401574 | 39 | 404593 | 38 |
| 401575 | 39 | 404594 | 38 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FIST

RESOURCES

RESOURCES

- Patient Ostomy Support & Resources
- Quick Links Guide
- Medicare Utilization Chart
- **Product Code Index**
- References

Request Product Sample

| SKU # | PG | SKU # | PG |
|-------------|----|-------------|----|
| 405181..... | 50 | 411657..... | 49 |
| 405446..... | 52 | 411658..... | 49 |
| 405447..... | 52 | 411659..... | 49 |
| 405448..... | 52 | 411800..... | 37 |
| 405449..... | 52 | 411801..... | 41 |
| 405450..... | 52 | 411802..... | 37 |
| 405451..... | 52 | 411803..... | 41 |
| 405452..... | 52 | 411804..... | 37 |
| 405453..... | 52 | 411805..... | 41 |
| 405456..... | 49 | 411806..... | 37 |
| 405457..... | 49 | 411807..... | 41 |
| 405458..... | 49 | 413153..... | 35 |
| 405467..... | 49 | 413154..... | 35 |
| 405468..... | 49 | 413155..... | 35 |
| 405473..... | 50 | 413156..... | 35 |
| 405474..... | 50 | 413157..... | 35 |
| 405475..... | 50 | 413159..... | 35 |
| 405476..... | 50 | 413160..... | 35 |
| 405477..... | 50 | 413161..... | 35 |
| 405478..... | 50 | 413162..... | 35 |
| 405479..... | 50 | 413163..... | 35 |
| 405480..... | 50 | 413164..... | 35 |
| 405481..... | 50 | 413165..... | 35 |
| 405482..... | 50 | 413166..... | 35 |
| 409268..... | 50 | 413167..... | 35 |
| 409269..... | 50 | 413168..... | 35 |
| 409270..... | 50 | 413177..... | 38 |
| 411264..... | 44 | 413178..... | 38 |
| 411265..... | 44 | 413179..... | 38 |
| 411266..... | 44 | 413180..... | 38 |
| 411267..... | 44 | 413181..... | 38 |
| 411288..... | 44 | 413182..... | 38 |
| 411289..... | 44 | 413183..... | 38 |
| 411290..... | 44 | 413184..... | 38 |
| 411291..... | 44 | 413185..... | 38 |
| 411292..... | 44 | 413186..... | 38 |
| 411489..... | 44 | 413187..... | 38 |
| 411490..... | 44 | 413188..... | 38 |
| 411491..... | 44 | 413417..... | 37 |
| 411492..... | 44 | 413418..... | 37 |
| 411493..... | 44 | 413419..... | 37 |
| 411631..... | 83 | 413420..... | 37 |
| 411632..... | 83 | 413421..... | 41 |
| 411633..... | 84 | 413422..... | 41 |
| 411637..... | 85 | 413423..... | 41 |
| 411638..... | 85 | 413424..... | 41 |
| 411639..... | 85 | 413437..... | 47 |
| 411640..... | 85 | 413438..... | 47 |
| 411641..... | 85 | 413439..... | 47 |
| 411642..... | 85 | 413509..... | 68 |

| SKU # | PG | SKU # | PG |
|-------------|----|-------------|----|
| 413510..... | 68 | 416740..... | 63 |
| 413515..... | 64 | 416741..... | 63 |
| 413516..... | 64 | 416742..... | 63 |
| 413517..... | 64 | 416743..... | 63 |
| 413518..... | 64 | 416744..... | 63 |
| 413519..... | 64 | 416745..... | 63 |
| 413520..... | 64 | 416746..... | 63 |
| 413521..... | 64 | 416747..... | 63 |
| 413522..... | 64 | 416748..... | 63 |
| 413523..... | 64 | 416749..... | 63 |
| 416403..... | 45 | 416750..... | 63 |
| 416406..... | 45 | 416751..... | 63 |
| 416409..... | 45 | 416774..... | 51 |
| 416412..... | 45 | 416775..... | 51 |
| 416414..... | 42 | 416776..... | 52 |
| 416415..... | 42 | 416777..... | 52 |
| 416416..... | 42 | 416778..... | 51 |
| 416417..... | 42 | 416779..... | 51 |
| 416418..... | 42 | 416781..... | 52 |
| 416419..... | 42 | 416782..... | 51 |
| 416420..... | 42 | 416783..... | 51 |
| 416421..... | 42 | 416795..... | 51 |
| 416422..... | 42 | 416796..... | 51 |
| 416423..... | 42 | 416798..... | 51 |
| 416424..... | 42 | 416799..... | 51 |
| 416472..... | 42 | 416901..... | 51 |
| 416700..... | 67 | 416902..... | 51 |
| 416701..... | 67 | 416903..... | 51 |
| 416703..... | 68 | 416908..... | 62 |
| 416704..... | 68 | 416975..... | 61 |
| 416706..... | 68 | 416976..... | 61 |
| 416707..... | 68 | 420695..... | 44 |
| 416709..... | 68 | 420696..... | 44 |
| 416710..... | 68 | 420697..... | 44 |
| 416712..... | 68 | 421025..... | 61 |
| 416713..... | 68 | 421026..... | 61 |
| 416715..... | 68 | 421027..... | 61 |
| 416718..... | 61 | 421028..... | 61 |
| 416719..... | 61 | 421029..... | 61 |
| 416721..... | 61 | 421033..... | 41 |
| 416724..... | 63 | 421034..... | 41 |
| 416725..... | 63 | 421035..... | 41 |
| 416728..... | 63 | 421039..... | 36 |
| 416729..... | 63 | 421040..... | 36 |
| 416732..... | 63 | 421041..... | 36 |
| 416733..... | 63 | 421074..... | 47 |
| 416736..... | 63 | 421454..... | 36 |
| 416737..... | 63 | 421456..... | 40 |
| 416738..... | 63 | 421458..... | 36 |
| 416739..... | 63 | 421460..... | 40 |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FIST

RESOURCES

RESOURCES

- Patient Ostomy Support & Resources
- Quick Links Guide
- Medicare Utilization Chart
- **Product Code Index**
- References

Request Product Sample

| SKU # | PG | SKU # | PG |
|--------|----|--------|-----|
| 421462 | 36 | 421853 | 61 |
| 421464 | 40 | 421854 | 63 |
| 421615 | 65 | 421855 | 63 |
| 421616 | 65 | 421856 | 63 |
| 421617 | 65 | 421857 | 63 |
| 421618 | 65 | 421858 | 61 |
| 421619 | 65 | 421859 | 63 |
| 421620 | 65 | 421860 | 63 |
| 421621 | 65 | 421861 | 63 |
| 421622 | 65 | 421862 | 63 |
| 421623 | 69 | 421893 | 45 |
| 421624 | 69 | 422159 | 65 |
| 421625 | 69 | 422163 | 97 |
| 421626 | 69 | 422363 | 66 |
| 421627 | 69 | 422364 | 66 |
| 421628 | 69 | 422365 | 66 |
| 421629 | 69 | 422539 | 66 |
| 421631 | 69 | 422540 | 66 |
| 421638 | 37 | 422541 | 66 |
| 421639 | 37 | 422545 | 66 |
| 421640 | 37 | 422546 | 66 |
| 421641 | 37 | 422551 | 70 |
| 421642 | 37 | 422552 | 70 |
| 421643 | 37 | 422553 | 70 |
| 421676 | 45 | 422557 | 70 |
| 421677 | 45 | 422558 | 70 |
| 421678 | 45 | 422559 | 70 |
| 421679 | 45 | 422560 | 70 |
| 421680 | 45 | 422562 | 70 |
| 421681 | 45 | 422564 | 70 |
| 421682 | 45 | 422565 | 70 |
| 421683 | 45 | 423283 | 92 |
| 421687 | 68 | 423284 | 92 |
| 421688 | 68 | 423285 | 92 |
| 421689 | 68 | 423286 | 91 |
| 421738 | 42 | 423287 | 92 |
| 421739 | 42 | 423288 | 91 |
| 421740 | 42 | 423289 | 100 |
| 421746 | 42 | 423290 | 100 |
| 421747 | 42 | 423391 | 99 |
| 421748 | 42 | 423392 | 91 |
| 421749 | 42 | 423548 | 78 |
| 421796 | 46 | 423549 | 78 |
| 421797 | 46 | 423550 | 79 |
| 421798 | 45 | 423551 | 79 |
| 421816 | 67 | 423552 | 79 |
| 421821 | 67 | 423553 | 79 |
| 421822 | 68 | 423554 | 77 |
| 421823 | 68 | 423555 | 77 |
| 421848 | 62 | 423556 | 77 |

| SKU # | PG | SKU # | PG |
|--------|-----|--------|-----|
| 423557 | 77 | 650830 | 73 |
| 423574 | 97 | 650831 | 73 |
| 423643 | 55 | 650832 | 73 |
| 423644 | 55 | 650833 | 73 |
| 423645 | 55 | 839001 | 93 |
| 423646 | 55 | 839002 | 93 |
| 423647 | 55 | 839003 | 95 |
| 423648 | 55 | 839004 | 95 |
| 423651 | 55 | 839005 | 93 |
| 423652 | 55 | 839006 | 93 |
| 423653 | 55 | 839010 | 94 |
| 423654 | 55 | 839250 | 105 |
| 423656 | 56 | 839251 | 105 |
| 423657 | 56 | 839252 | 105 |
| 423658 | 56 | 839253 | 106 |
| 423659 | 56 | 839261 | 105 |
| 423660 | 56 | 839262 | 105 |
| 423661 | 56 | 839264 | 106 |
| 423664 | 56 | 839265 | 106 |
| 423665 | 56 | 839266 | 106 |
| 423666 | 56 | 839267 | 105 |
| 423667 | 56 | 839268 | 106 |
| 423676 | 57 | 839280 | 106 |
| 423677 | 57 | TR105 | 98 |
| 423678 | 57 | | |
| 423689 | 57 | | |
| 423690 | 57 | | |
| 423691 | 57 | | |
| 423704 | 58 | | |
| 423705 | 58 | | |
| 423712 | 58 | | |
| 423713 | 58 | | |
| 423714 | 58 | | |
| 423715 | 58 | | |
| 423716 | 58 | | |
| 423717 | 58 | | |
| 423718 | 58 | | |
| 423722 | 59 | | |
| 423723 | 59 | | |
| 423725 | 59 | | |
| 423726 | 59 | | |
| 423727 | 59 | | |
| 423779 | 92 | | |
| 423780 | 100 | | |
| 423783 | 100 | | |
| 423784 | 92 | | |
| 432558 | 77 | | |
| 432559 | 77 | | |
| 650828 | 73 | | |
| 650829 | 73 | | |

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FIST

RESOURCES

RESOURCES

- Patient Ostomy Support & Resources
- Quick Links Guide
- Medicare Utilization Chart
- **Product Code Index**
- References

Request Product Sample


1. Bulkeley, J.E., McMullen, C.K., Grant, M., Wendel, C., Hornbrook, M.C., Krouse, R.S. (2018). Ongoing Ostomy Self-Care Challenges of Long-Term Rectal Cancer Survivors. Supportive Care in Cancer, 26 (11), 3933-3939.
2. Consumer survey following me+ nurse consultation, N=166, Convatec, Data on File.
3. Consumer survey following telehealth nurse consultation, Convatec & Corstrata, 2023, data on file.
4. Hoeflok J, et al. A prospective multicenter evaluation of a moldable stoma skin barrier. Ostomy Wound Manage 2009; 55: 62-69.
5. Chaumier D, et al. An evaluation of the peristomal skin condition in ostomates using moldable skin barriers. 44th Annual WOCN conference, Charlotte, NC, USA; June 11, 2012.
6. Szewczyk MT, et al. The effects of using a moldable skin barrier on peristomal skin condition in persons with an ostomy: results of a prospective, observational, multinational study. Ostomy Wound Manage 2014; 60: 16-26.
7. Ratliff, C., Goldbert, M., Jaszarowski, K, McNichol, L., Pittman, J., Gray, M. (2021). Peristomal Skin Health, A WOCN Society Consensus Conference. Journal of Wound, Ostomy Continence Nursing, 2021;48(3):219-231.
8. Maria Teresa Szewczyk, MD, PhD; Grazyna Majewska, RN, ETN; Mary V. Cabral, MS, FNP-BC, CWOCN-AP; and Karin Hölzel-Piontek, RN; The Effects of Using a Moldable Skin Barrier on Peristomal Skin Condition in Persons with an Ostomy: Results of a Prospective, Observational, Multinational Study, Ostomy Wound Management 2014;60(12):16-26.
9. Clinical Evaluation of Sur Fit Natura Cut to Fit Skin Wafers Gillian Seaton. Technical data on file, reference document WHRI3940 TA238 2014.
10. Rolstad, BS. Principles and techniques in the use of convexity. Ostomy Wound Manage. 1996. 42(1); 24-32.
11. Pontieri, Lewis V. Basics of ostomy care. Medsurg Nurs. 2006; 15(4): 199-202.
12. Hanley, J. Convex stoma appliances: are we getting it right? British Journal of Nursing (Stoma Supplement). 2013; 22(16).
13. Turnbull, GB. The convexity controversy. Ostomy Wound Manage. 2003; 49(1): 16- www.o-wm.com/article/1203. Accessed 29 Feb, 2016.
14. DRG data 2015, in T5 markets: Japan, UK, USA, Germany, China.
15. Boyd, K. Clinical Protocols for Stoma Care - Use of convex appliances. Nursing Standard. 2004; 18(20): 37-38.
16. Burch, 2011: Management of Stoma Complications. NursingTimes.net. 10 Nov 2011.
17. McNichol L, Cobb T, Depaive Y, Quigley M, Smitka K, Gray, M. Characteristics of Convex Skin Barriers and Clinical Application: Results of an International Consensus Panel. J Wound Ostomy Continence Nurs. 2021;48(6):524-532
18. Voegeli D. LBF® "No-Sting" Barrier Wipes: Skin Care Using Advanced Silicone Technology. British Journal of Nursing, 2008; Vol. 17: No. 7
19. Data on file. Convatec Inc.

Disclaimer:

The reimbursement information provided by Convatec is intended to provide general information relevant to coding and reimbursement of Convatec's products only. Coverage and payment policies for the same insurer can vary from one region to another and may change from time to time because of ongoing changes in government and insurance industry rules and regulations. Therefore, please confirm HCPCS codes with your local DME-MAC, private insurer, or Medicaid agency before processing claims. Convatec does not guarantee coverage or payment of its products listed herein.



Leakage problem? Problem solved.


[LEARN MORE](#)

visit [convatec.com](https://www.convatec.com)
call 877-585-0470 or email us at
concierge@Convatec.com

Introducing EsteemBody™ with Leak Defense™

Leak Defense™ combines gold-standard adhesives with a comprehensive soft convexity range designed to adapt to the body for a secure, longer-lasting seal.

- ✓ Our modern 8-shaped pouch is designed to sag less as it fills.*
- ✓ As pioneers in hydrocolloids, Convatec's adhesives were created to help users get it right the first time every time.
- ✓ Our range of soft convexity combinations were designed to adapt to the needs of unique body and stoma types.

TABLE OF CONTENTS

[me+™ SUPPORT](#)
[EDUCATION](#)
[TWO-PIECE SYSTEMS](#)
[ONE-PIECE SYSTEMS](#)
[ACUTE KITS](#)
[PEDIATRIC PRODUCTS](#)
[ACCESSORIES](#)
[RESOURCES](#)

- Patient Ostomy Support & Resources
- Quick Links Guide
- Medicare Utilization Chart
- Product Code Index
- **References**


[WOUND & FIST](#)
[RESOURCES](#)
[Request Product Sample](#)

me+™ is here for you and your patients

With our ostomy care apps, you
and your patients can access
support wherever you are using
your smart phone or tablet.




For general product information,
call 1-800-422-8811
or visit Convatec.com/ostomy



LEARN MORE

about the
Convatec
Ostomy Nurse
Solutions™ App



LEARN MORE

about the
My Ostomy
Journey™ App

TABLE OF CONTENTS

me+™ SUPPORT

EDUCATION

TWO-PIECE SYSTEMS

ONE-PIECE SYSTEMS

ACUTE KITS

PEDIATRIC PRODUCTS

ACCESSORIES

WOUND & FISTULA POUCHES

RESOURCES

Request Product Sample